

Package: raveio (via r-universe)

October 20, 2026

Type Package

Title File-System Toolbox for RAVE Project

Version 0.9.0.101

Language en-US

Description Includes multiple cross-platform read/write interfaces for 'RAVE' project. 'RAVE' stands for ``R analysis and visualization of human intracranial electroencephalography data". The whole project aims at providing powerful free-source package that analyze brain recordings from patients with electrodes placed on the cortical surface or inserted into the brain. 'raveio' as part of this project provides tools to read/write neurophysiology data from/to 'RAVE' file structure, as well as several popular formats including 'EDF(+)', 'Matlab', 'BIDS-iEEG', and 'HDF5', etc. Documentation and examples about 'RAVE' project are provided at <<https://rave.wiki>>, and the paper by John F. Magnotti, Zhengjia Wang, Michael S. Beauchamp (2020) <[doi:10.1016/j.neuroimage.2020.117341](https://doi.org/10.1016/j.neuroimage.2020.117341)>. For applications such as electrode localization, please check Zhengjia Wang, John F. Magnotti, Xiang Zhang, Michael S. Beauchamp (2023) <[doi:10.1523/ENEURO.0328-23.2023](https://doi.org/10.1523/ENEURO.0328-23.2023)>, or see 'citation(` `raveio)'' for details.

BugReports <https://github.com/beauchamplab/raveio/issues>

URL <https://beauchamplab.github.io/raveio/>

License MIT + file LICENSE

Encoding UTF-8

RoxygenNote 7.3.3

SystemRequirements little-endian platform

Imports utils, stats, data.table, dplyr, edfReader (>= 1.2.1), filearray (>= 0.1.3), ieegio, R6, ravepipeline, ravetools, rpyANTs, targets (>= 0.8.0), threeBrain (>= 0.2.5)

Suggests callr (>= 3.7.0), stringr (>= 1.4.0), globals, jsonlite, testthat, shiny, rpy2, freesurferformats, readNSx, gsignal

Config/pak/sysreqs cmake libfftw3-dev git libglpk-dev make libicu-dev libpng-dev libuv1-dev libxml2-dev libssl-dev python3 libzmq3-dev zlib1g-dev

Repository <https://dipterix.r-universe.dev>

Date/Publication 2026-04-28 19:40:51 UTC

RemoteUrl <https://github.com/beauchamplab/raveio>

RemoteRef HEAD

RemoteSha efd87d4ba61e9ff5fbb8ee19f9a3bb8d00591b87

Contents

ants_coreg	4
ants_preprocessing	6
archive_subject	6
as_rave_project	8
as_rave_subject	9
as_rave_unit	9
as_yael_process	10
Auxiliary_electrode	11
backup_file	13
BlackrockFile	14
cache_path	16
cache_to_filearray	17
catgl	19
cmd_run_3dAllineate	20
cmd_run_yael_preprocess	23
collapse_power	26
collapse2	28
compose_channel	29
convert_blackrock	31
convert_electrode_table_to_bids	32
dir_create2	33
ECoGTensor	33
export_table	35
find_path	36
generate_reference	37
get_projects	38
get_val2	39
h5_names	40
h5_valid	40
import_electrode_table	41
ingest_regressor	42
install_subject	43
is.blank	44
is.zerolenth	44
is_on_cran	45
is_valid_ish	45

join_tensors	47
lapply_async	48
LFP_electrode	49
LFP_reference	54
load_bids_ieeg_header	58
load_fst_or_h5	59
load_h5	60
load_meta2	61
load_yaml	62
mgh_to_nii	62
new_electrode	63
niftyreg_coreg	64
power_baseline	66
prepare_subject_bare0	68
progress_with_logger	72
py_nipy_coreg	73
rave-raw-validation	74
rave-server	76
rave_brain	77
rave_command_line_path	78
rave_directories	79
rave_export	80
rave_import	81
rave_subject_format_conversion	83
RAVEAbstarctElectrode	84
RAVEEpoch	87
raveio-constants	90
RAVEMetaSubject	91
RAVEPreprocessSettings	93
RAVEProject	96
RAVESubject	98
read-brainvision-eeg	103
read-write-fst	105
read_csv_ieeg	105
read_edf_header	106
read_edf_signal	107
read_mat	107
read_nsx_nev	109
safe_read_csv	110
safe_write_csv	111
save_h5	112
save_json	113
save_meta2	114
save_yaml	115
Spike_electrode	116
Tensor	118
test_hdspeed	123
time_diff2	124

transform_point_to_template	125
url_neurosynth	126
validate_subject	127
validate_time_window	128
voltage_baseline	129
with_future_parallel	132
YAELProcess	133

Index 140

ants_coreg	<i>Register 'CT' or 'MR' images via 'ANTs'</i>
------------	--

Description

ants_coreg aligns 'CT' to 'MR' images; ants_mri_to_template aligns native 'MR' images to group templates

Usage

```
ants_coreg(
  ct_path,
  mri_path,
  coreg_path = NULL,
  reg_type = c("DenseRigid", "Rigid", "SyN", "Affine", "TRSAA", "SyNCC", "SyNOnly"),
  aff_metric = c("mattes", "meansquares", "GC"),
  syn_metric = c("mattes", "meansquares", "demons", "CC"),
  verbose = TRUE,
  ...
)

cmd_run_ants_coreg(
  subject,
  ct_path,
  mri_path,
  reg_type = c("DenseRigid", "Rigid", "SyN", "Affine", "TRSAA", "SyNCC", "SyNOnly"),
  aff_metric = c("mattes", "meansquares", "GC"),
  syn_metric = c("mattes", "meansquares", "demons", "CC"),
  verbose = TRUE,
  dry_run = FALSE
)

ants_mri_to_template(
  subject,
  template_subject = getOption("threeBrain.template_subject", "N27"),
  preview = FALSE,
  verbose = TRUE,
  ...
)
```

```

)

cmd_run_ants_mri_to_template(
    subject,
    template_subject = getOption("threeBrain.template_subject", "N27"),
    verbose = TRUE,
    dry_run = FALSE
)

ants_morph_electrode(subject, preview = FALSE, dry_run = FALSE)

```

Arguments

ct_path, mri_path	absolute paths to 'CT' and 'MR' image files
coreg_path	registration path, where to save results; default is the parent folder of ct_path
reg_type	registration type, choices are 'DenseRigid', 'Rigid', 'Affine', 'SyN', 'TRSAA', 'SyNCC', 'SyNOnly', or other types; see ants_registration
aff_metric	cost function to use for linear or 'affine' transform
syn_metric	cost function to use for 'SyN' transform
verbose	whether to verbose command; default is true
...	other arguments passed to ants_registration
subject	'RAVE' subject
dry_run	whether to dry-run the script and to print out the command instead of executing the code; default is false
template_subject	template to map 'MR' images
preview	whether to preview results; default is false

Value

Aligned 'CT' will be generated at the coreg_path path:

- 'ct_in_t1.nii.gz' aligned 'CT' image; the image is also re-sampled into 'MRI' space
- 'transform.yaml' transform settings and outputs
- 'CT_IJK_to_MR_RAS.txt' transform matrix from volume 'IJK' space in the original 'CT' to the 'RAS' anatomical coordinate in 'MR' scanner; 'affine' transforms only
- 'CT_RAS_to_MR_RAS.txt' transform matrix from scanner 'RAS' space in the original 'CT' to 'RAS' in 'MR' scanner space; 'affine' transforms only

ants_preprocessing	<i>Process 'T1' weighted 'MRI' using 'ANTs'</i>
--------------------	---

Description

This function is soft-deprecated. Use [yael_preprocess](#) instead.

Usage

```
ants_preprocessing(
  work_path,
  image_path,
  resample = FALSE,
  verbose = TRUE,
  template_subject = raveio_getopt("threeBrain_template_subject")
)
```

Arguments

work_path	working directory, all intermediate images will be stored here
image_path	input image path
resample	whether to resample the input image before processing
verbose	whether to verbose the processing details
template_subject	template mapping, default is derived from raveio_getopt

Value

Nothing. All images are saved to work_path

archive_subject	<i>Archive and share a subject</i>
-----------------	------------------------------------

Description

Archive and share a subject

Usage

```
archive_subject(
  subject,
  path,
  includes = c("original_signals", "processed_data", "rave_imaging", "pipelines", "notes",
    "user_generated"),
  config = list(),
  work_path = NULL,
  zip_flags = NULL
)
```

Arguments

subject	'RAVE' subject to archive
path	path to a zip file to store; if missing or empty, then the path will be automatically created
includes	data to include in the archive; default includes all (original raw signals, processed signals, imaging files, stored pipelines, notes, and user-generated exports)
config	a list of configurations, including changing subject code, project name, or to exclude cache data; see examples
work_path	temporary working path where files are copied; default is temporary path. Set this variable explicitly when temporary path is on external drives (for example, users have limited storage on local drives and cannot hold the entire subject)
zip_flags	zip flags

Examples

```
# This example requires you to install demo subject

## Not run:

# Basic usage
path <- archive_subject('demo/DemoSubject')

# clean up
unlink(path)

# Advanced usage: include all the original signals
# and processed data, no cache data, re-name to
# demo/DemoSubjectLite
path <- archive_subject(
  'demo/DemoSubject',
  includes = c("original_signals", "processed_data"),
  config = list(
    rename = list(
      project_name = "demo",
      subject_code = "DemoSubjectLite"
```

```

    ),
    original_signals = list(
        # include all raw signals
        include_all = TRUE
    ),
    processed_data = list(
        include_cache = FALSE
    )
)
)

# Clean up temporary zip file
unlink(path)

## End(Not run)

```

as_rave_project	Convert character to RAVEProject instance
-----------------	---

Description

Convert character to [RAVEProject](#) instance

Usage

```
as_rave_project(project, ...)
```

Arguments

project	character project name
...	passed to other methods

Value

A [RAVEProject](#) instance

See Also

[RAVEProject](#)

as_rave_subject	Get RAVESubject instance from character
-----------------	---

Description

Get [RAVESubject](#) instance from character

Usage

```
as_rave_subject(subject_id, strict = TRUE, reload = TRUE)
```

Arguments

subject_id	character in format "project/subject"
strict	whether to check if subject directories exist or not
reload	whether to reload (update) subject information, default is true

Value

[RAVESubject](#) instance

See Also

[RAVESubject](#)

as_rave_unit	Convert numeric number into print-friendly format
--------------	---

Description

Convert numeric number into print-friendly format

Usage

```
as_rave_unit(x, unit, label = "")
```

Arguments

x	numeric or numeric vector
unit	the unit of x
label	prefix when printing x

Value

Still numeric, but print-friendly class

Examples

```

sp <- as_rave_unit(1024, 'GB', 'Hard-disk space is ')
print(sp, digits = 0)

sp - 12

as.character(sp)

as.numeric(sp)

# Vectorize
sp <- as_rave_unit(c(500,200), 'MB/s', c('Writing: ', 'Reading: '))
print(sp, digits = 0, collapse = '\n')

```

as_yael_process	Create a 'Yael' imaging processing instance
-----------------	---

Description

Image registration across different modals. Normalize brain 'T1'-weighted 'MRI' to template brain and generate subject-level atlas files.

Usage

```
as_yael_process(subject)
```

Arguments

subject character (subject code, or project name with subject code), or [RAVESubject](#) instance.

Value

A processing instance, see [YaelProcess](#)

Examples

```

library(raveio)
process <- as_yael_process("testtest2")

# This example requires extra demo data & settings.
## Not run:

# Import and set original T1w MRI and CT
process$set_input_image("/path/to/T1w_MRI.nii", type = "T1w")
process$set_input_image("/path/to/CT.nii.gz", type = "CT")

# Co-register CT to MRI

```

```

process$register_to_T1w(image_type = "CT")

# Morph T1w MRI to 0.5 mm^3 MNI152 template
process$map_to_template(
  template_name = "mni_icbm152_nlin_asym_09b",
  native_type = "T1w"
)

## End(Not run)

```

Auxiliary_electrode *Class definition for auxiliary channels*

Description

Class definition for auxiliary channels

Class definition for auxiliary channels

Value

If simplify is enabled, and only one block is loaded, then the result will be a vector (type="voltage") or a matrix (others), otherwise the result will be a named list where the names are the blocks.

Super class

`raveio::RAVEAbstarctElectrode` -> Auxiliary_electrode

Active bindings

`h5_fname` 'HDF5' file name

`valid` whether current electrode is valid: subject exists and contains current electrode or reference;
subject electrode type matches with current electrode type

`raw_sample_rate` voltage sample rate

`preprocess_info` preprocess information

`voltage_file` path to voltage 'HDF5' file

Methods

Public methods:

- `Auxiliary_electrode$print()`
- `Auxiliary_electrode$set_reference()`
- `Auxiliary_electrode$new()`
- `Auxiliary_electrode$.load_noref_voltage()`

- `Auxiliary_electrode$.load_raw_voltage()`
- `Auxiliary_electrode$load_data()`
- `Auxiliary_electrode$load_blocks()`
- `Auxiliary_electrode$clear_cache()`
- `Auxiliary_electrode$clear_memory()`
- `Auxiliary_electrode$clone()`

Method `print()`: print electrode summary

Usage:

`Auxiliary_electrode$print()`

Method `set_reference()`: set reference for current electrode

Usage:

`Auxiliary_electrode$set_reference(reference)`

Arguments:

reference either NULL or LFP_electrode instance

Method `new()`: constructor

Usage:

`Auxiliary_electrode$new(subject, number, quiet = FALSE)`

Arguments:

subject, number, quiet see constructor in [RAVEAbstarctElectrode](#)

Method `.load_noref_voltage()`: load non-referenced voltage (internally used)

Usage:

`Auxiliary_electrode$.load_noref_voltage(reload = FALSE)`

Arguments:

reload whether to reload cache

srate voltage signal sample rate

Method `.load_raw_voltage()`: load raw voltage (no process)

Usage:

`Auxiliary_electrode$.load_raw_voltage(reload = FALSE)`

Arguments:

reload whether to reload cache

Method `load_data()`: method to load electrode data

Usage:

`Auxiliary_electrode$load_data(type = c("raw-voltage", "voltage"))`

Arguments:

type data type such as "power", "phase", "voltage", "wavelet-coefficient", and "raw-voltage".

For "power", "phase", and "wavelet-coefficient", 'Wavelet' transforms are required.

For "voltage", 'Notch' filters must be applied. All these types except for "raw-voltage" will be referenced. For "raw-voltage", no reference will be performed since the data will be the "raw" signal (no processing).

Method `load_blocks()`: load electrode block-wise data (with no reference), useful when epoch is absent

Usage:

```
Auxiliary_electrode$load_blocks(
  blocks,
  type = c("raw-voltage", "voltage"),
  simplify = TRUE
)
```

Arguments:

`blocks` session blocks

`type` data type such as "power", "phase", "voltage", "raw-voltage" (with no filters applied, as-is from imported), "wavelet-coefficient". Note that if type is "raw-voltage", then the data only needs to be imported; for "voltage" data, 'Notch' filters must be applied; for all other types, 'Wavelet' transforms are required.

`simplify` whether to simplify the result

Method `clear_cache()`: method to clear cache on hard drive

Usage:

```
Auxiliary_electrode$clear_cache(...)
```

Arguments:

`...` ignored

Method `clear_memory()`: method to clear memory

Usage:

```
Auxiliary_electrode$clear_memory(...)
```

Arguments:

`...` ignored

Method `clone()`: The objects of this class are cloneable with this method.

Usage:

```
Auxiliary_electrode$clone(deep = FALSE)
```

Arguments:

`deep` Whether to make a deep clone.

backup_file

Back up and rename the file or directory

Description

Back up and rename the file or directory

Usage

```
backup_file(path, remove = FALSE, quiet = FALSE)
```

Arguments

path	path to a file or a directory
remove	whether to remove the original path; default is false
quiet	whether not to verbose the messages; default is false

Value

FALSE if nothing to back up, or the back-up path if path exists

Examples

```
path <- tempfile()
file.create(path)

path2 <- backup_file(path, remove = TRUE)

file.exists(c(path, path2))
unlink(path2)
```

BlackrockFile

Class definition to load data from 'BlackRock' 'Micro-systems' files

Description

Currently only supports minimum file specification version 2.3. Please contact the package maintainer or 'RAVE' team if older specifications are needed

Value

- absolute file path
- absolute file paths
- nothing
- a data frame
- a list of spike 'waveform' (without normalization)
- a normalized numeric vector (analog signals with 'uV' as the unit)

Public fields

- block character, session block ID

Active bindings

base_path absolute base path to the file
 version 'NEV' specification version
 electrode_table electrode table
 sample_rate_nev_timestamp sample rate of 'NEV' data packet time-stamps
 has_nsx named vector of 'NSx' availability
 recording_duration recording duration of each 'NSx'
 sample_rates sampling frequencies of each 'NSx' file

Methods**Public methods:**

- `BlackrockFile#print()`
- `BlackrockFile$new()`
- `BlackrockFile$nev_path()`
- `BlackrockFile$nsx_paths()`
- `BlackrockFile$refresh_data()`
- `BlackrockFile$get_epoch()`
- `BlackrockFile$get_waveform()`
- `BlackrockFile$get_electrode()`
- `BlackrockFile$clone()`

Method `print()`: print user-friendly messages

Usage:

`BlackrockFile#print()`

Method `new()`: constructor

Usage:

`BlackrockFile$new(path, block, nev_data = TRUE)`

Arguments:

path the path to 'BlackRock' file, can be with or without file extensions

block session block ID; default is the file name

nev_data whether to load comments and 'waveforms'

Method `nev_path()`: get 'NEV' file path

Usage:

`BlackrockFile$nev_path()`

Method `nsx_paths()`: get 'NSx' file paths

Usage:

`BlackrockFile$nsx_paths(which = NULL)`

Arguments:

which which signal file to get, or NULL to return all available paths, default is NULL; must be integers

Method refresh_data(): refresh and load 'NSx' data

Usage:

BlackrockFile\$refresh_data(force = FALSE, verbose = TRUE, nev_data = FALSE)

Arguments:

force whether to force reload data even if the data has been loaded and cached before

verbose whether to print out messages when loading

nev_data whether to refresh 'NEV' extended data; default is false

Method get_epoch(): get epoch table from the 'NEV' comment data packet

Usage:

BlackrockFile\$get_epoch()

Method get_waveform(): get 'waveform' of the spike data

Usage:

BlackrockFile\$get_waveform()

Method get_electrode(): get electrode data

Usage:

BlackrockFile\$get_electrode(electrode, nstype = NULL)

Arguments:

electrode integer, must be a length of one

nstype which signal bank, for example, 'ns3', 'ns5'

Method clone(): The objects of this class are cloneable with this method.

Usage:

BlackrockFile\$clone(deep = FALSE)

Arguments:

deep Whether to make a deep clone.

cache_path

Manipulate cached data on the file systems

Description

Manipulate cached data on the file systems

Usage

cache_root(check = FALSE)

clear_cached_files(subject_code, quiet = FALSE)

Arguments

check	whether to ensure the cache root path
subject_code	subject code to remove; default is missing. If subject_code is provided, then only this subject-related cache files will be removed.
quiet	whether to suppress the message

Details

'RAVE' intensively uses cache files. If running on personal computers, the disk space might be filled up very quickly. These cache files are safe to remove if there is no 'RAVE' instance running. Function `clear_cached_files` is designed to remove these cache files. To run this function, please make sure that all 'RAVE' instances are shutdown.

Value

`cache_root` returns the root path that stores the 'RAVE' cache data; `clear_cached_files` returns nothing

Examples

```
cache_root()
```

cache_to_filearray	<i>Generate and automatically cache a file array</i>
--------------------	--

Description

Avoid repeating yourself

Usage

```
cache_to_filearray(
  fun,
  filebase,
  globals,
  dimension,
  type = "auto",
  partition_size = 1L,
  verbose = FALSE,
  ...
)
```

Arguments

fun	function that can be called with no mandatory arguments; the result should be in a matrix or an array
filebase	where to store the array
globals	names of variables such that any changes should trigger a new evaluation of fun. This argument is highly recommended to be set explicitly (with atomic variables) though the function automatically calculates the global variables
dimension	what is the supposed dimension, default is automatically calculated from array. If specified explicitly and the file array dimension is inconsistent, a new calculation will be triggered.
type	file array type, default is "auto"; can be explicitly specified; see filearray_create . Inconsistent type will trigger a new calculation.
partition_size	file array partition size; default is 1; set it to NA to generate it automatically. Notice inconsistent partition size will not trigger calculation if the key variables remain the same
verbose	whether to verbose debug information
...	passed to findGlobals

Examples

```

c <- 2
b <- list(d = matrix(1:9,3))
filebase <- tempfile()

f <- function() {
  message("New calculation")
  re <- c + b$d
  dimnames(re) <- list(A=1:3, B = 11:13)

  # `extra` attribute will be saved
  attr(re, "extra") <- "extra meta data"
  re
}

# first time running
arr <- cache_to_filearray( f, filebase = filebase )

# cached, no re-run
arr <- cache_to_filearray( f, filebase = filebase )

# file array object
arr

# read into memory
arr[]

# read extra data

```

```

arr$get_header("extra")

# get digest results
arr$get_header("raveio::filearray_cache")

## Clean up this example
unlink(filebase, recursive = TRUE)

```

catgl	<i>Print colored messages</i>
-------	-------------------------------

Description

Print colored messages

Usage

```
catgl(..., .envir = parent.frame(), level = "DEBUG", .pal, .capture = FALSE)
```

Arguments

..., .envir	passed to glue
level	passed to cat2
.pal	see pal in cat2
.capture	logical, whether to capture message and return it without printing

Details

The level has order that sorted from low to high: "DEBUG", "DEFAULT", "INFO", "WARNING", "ERROR", "FATAL". Each different level will display different colors and icons before the message. You can suppress messages with certain levels by setting 'raveio' options via `raveio_setopt('verbose_level', <level>)`. Messages with levels lower than the threshold will be muffled. See examples.

Value

The message as characters

Examples

```

# ----- Basic Styles -----

# Temporarily change verbose level for example
raveio_setopt('verbose_level', 'DEBUG', FALSE)

# debug
catgl('Debug message', level = 'DEBUG')

```

```

# default
catgl('Default message', level = 'DEFAULT')

# info
catgl('Info message', level = 'INFO')

# warning
catgl('Warning message', level = 'WARNING')

# error
catgl('Error message', level = 'ERROR')

try({
  # fatal, will call stop and raise error
  catgl('Error message', level = 'FATAL')
}, silent = TRUE)

# ----- Muffle messages -----

# Temporarily change verbose level to 'WARNING'
raveio_setopt('verbose_level', 'WARNING', FALSE)

# default, muffled
catgl('Default message')

# message printed for level >= Warning
catgl('Default message', level = 'WARNING')
catgl('Default message', level = 'ERROR')

```

cmd_run_3dAllineate *External shell commands for 'RAVE'*

Description

These shell commands are for importing 'DICOM' images to 'Nifti' format, reconstructing cortical surfaces, and align' the CT' to 'MRI'. The commands are only tested on 'MacOS' and 'Linux'. On 'Windows' machines, please use the 'WSL2' system.

Usage

```

cmd_run_3dAllineate(
  subject,
  mri_path,
  ct_path,
  overwrite = FALSE,
  command_path = NULL,

```

```
    dry_run = FALSE,
    verbose = dry_run
)

cmd_execute(
  script,
  script_path,
  command = "bash",
  dry_run = FALSE,
  backup = TRUE,
  args = NULL,
  ...
)

cmd_run_r(
  expr,
  quoted = FALSE,
  verbose = TRUE,
  dry_run = FALSE,
  log_file = tempfile(),
  script_path = tempfile(),
  ...
)

cmd_run_dcm2niix(
  subject,
  src_path,
  type = c("MRI", "CT"),
  merge = c("Auto", "No", "Yes"),
  float = c("Yes", "No"),
  crop = c("No", "Yes", "Ignore"),
  overwrite = FALSE,
  command_path = NULL,
  dry_run = FALSE,
  verbose = dry_run
)

cmd_run_flirt(
  subject,
  mri_path,
  ct_path,
  dof = 6,
  cost = c("mutualinfo", "leastsq", "normcorr", "corratio", "normmi", "labeldiff", "bbr"),
  search = 90,
  searchcost = c("mutualinfo", "leastsq", "normcorr", "corratio", "normmi", "labeldiff",
    "bbr"),
  overwrite = FALSE,
  command_path = NULL,
```

```

    dry_run = FALSE,
    verbose = dry_run
)

cmd_run_recon_all(
    subject,
    mri_path,
    args = c("-all", "-autorecon1", "-autorecon2", "-autorecon3", "-autorecon2-cp",
             "-autorecon2-wm", "-autorecon2-pial"),
    work_path = NULL,
    overwrite = FALSE,
    command_path = NULL,
    dry_run = FALSE,
    verbose = dry_run
)

cmd_run_recon_all_clinical(
    subject,
    mri_path,
    work_path = NULL,
    overwrite = FALSE,
    command_path = NULL,
    dry_run = FALSE,
    verbose = dry_run,
    ...
)

```

Arguments

subject	characters or a RAVESubject instance
mri_path	the absolute to 'MRI' volume; must in 'Nifti' format
ct_path	the absolute to 'CT' volume; must in 'Nifti' format
overwrite	whether to overwrite existing files; default is false
command_path	command line path if 'RAVE' cannot find the command binary files
dry_run	whether to run in dry-run mode; under such mode, the shell command will not execute. This is useful for debugging scripts; default is false
verbose	whether to print out the command script; default is true under dry-run mode, and false otherwise
script	the shell script
script_path	path to run the script
command	which command to invoke; default is 'bash'
backup	whether to back up the script file immediately; default is true
args	further arguments in the shell command, especially the 'FreeSurfer' reconstruction command
...	passed to system2 , or additional arguments

expr	expression to run as command
quoted	whether expr is quoted; default is false
log_file	where should log file be stored
src_path	source of the 'DICOM' or 'Nifti' image (absolute path)
type	type of the 'DICOM' or 'Nifti' image; choices are 'MRI' and 'CT'
merge, float, crop	'dcm2niix' conversion arguments; ignored when the source is in 'Nifti' format
dof, cost, search, searchcost	parameters used by 'FSL' 'flirt' command; see their documentation for details
work_path	work path for 'FreeSurfer' command;

Value

A list of data containing the script details:

script script details
script_path where the script should/will be saved
dry_run whether dry-run mode is turned on
log_file path to the log file
execute a function to execute the script

cmd_run_yael_preprocess

Process brain images for 'Yael'

Description

Aligns 'T1w' with other image types; normalizes 'T1w' 'MRI' to 'MNI152' templates via symmetric non-linear morphs. Create brain custom atlases from templates.

Usage

```
cmd_run_yael_preprocess(
  subject_code,
  t1w_path = NULL,
  ct_path = NULL,
  t2w_path = NULL,
  fgatir_path = NULL,
  preopct_path = NULL,
  flair_path = NULL,
  t1w_contrast_path = NULL,
  register_reversed = FALSE,
  normalize_template = "mni_icbm152_nlin_asym_09b",
```

```

    run_recon_all = TRUE,
    dry_run = FALSE,
    verbose = TRUE
)

yael_preprocess(
  subject_code,
  t1w_path = NULL,
  ct_path = NULL,
  t2w_path = NULL,
  fgatir_path = NULL,
  preopct_path = NULL,
  flair_path = NULL,
  t1w_contrast_path = NULL,
  register_policy = c("auto", "all"),
  register_reversed = FALSE,
  normalize_template = "mni_icbm152_nlin_asym_09a",
  normalize_policy = c("auto", "all"),
  normalize_back = ifelse(length(normalize_template) >= 1, normalize_template[[1]], NA),
  atlases = list(),
  add_surfaces = FALSE,
  verbose = TRUE
)

```

Arguments

subject_code	'RAVE' subject code
t1w_path	(required) 'T1' weighted 'MRI' path
ct_path	(optional in general but mandatory for electrode localization) post-surgery 'CT' path
t2w_path	(optional) 'T2' weighted 'MRI' path
fgatir_path	(optional) 'fGATIR' (fast gray-matter acquisition 'T1' inversion recovery) image path
preopct_path	(optional) pre-surgery 'CT' path
flair_path	(optional) 'FLAIR' (fluid-attenuated inversion recovery) image path
t1w_contrast_path	(optional) 'T1' weighted 'MRI' with contrast (usually used to show the blood vessels)
register_reversed	direction of the registration; FALSE (default) registers other images (such as post-surgery 'CT' to 'T1'); set to FALSE if you would like the 'T1' to be registered into other images. Since 'Yael' does not re-sample the images, there is no essential difference on the final registration results
normalize_template	names of the templates which the native 'T1' images will be normalized into
run_recon_all	whether to run 'FreeSurfer' reconstruction; default is true

dry_run	whether to dry-run the script and check if error exists before actually execute the scripts.
verbose	whether to print out the progress; default is TRUE
register_policy	whether images should be registered with 'T1w' image; default is "auto": automatically run registration algorithm if missing; alternative is "all": force the registration algorithm even if mapping files exist
normalize_policy	normalization policy; similar to register_policy but is applied to normalization. Default is "auto": automatically run normalization when the mapping is missing, and skip if exists; alternative is "all": force to run the normalization.
normalize_back	length of one (select from normalize_template), which template is to be used to generate native brain mask and transform matrices
atlases	a named list: the names must be template names from normalize_template and the values must be directories of atlases of the corresponding templates (see 'Examples').
add_surfaces	Whether to add surfaces for the subject; default is FALSE. The surfaces are created by reversing the normalization from template brain, hence the results will not be accurate. Enable this option only if cortical surface estimation is not critical.

Value

Nothing, a subject imaging folder will be created under 'RAVE' raw folder

Examples

```
## Not run:

# For T1 preprocessing only
yael_preprocess(
  subject_code = "patient01",
  t1w_path = "/path/to/T1.nii or T1.nii.gz",

  # normalize T1 to all 2009 MNI152-Asym brains (a,b,c)
  normalize_template = c(
    "mni_icbm152_nlin_asym_09b"
  ),

  # only normalize if not exists
  normalize_policy = "auto",

  # use MNI152b to create native processing folder
  normalize_back = "mni_icbm152_nlin_asym_09b",

  # Atlases generated from different templates have different
  # coordinates, hence both folder path and template names must be
  # provided
```

```

    atlases = list(
      mni_icbm152_nlin_asym_09b = "/path/to/atlas/folder1",
      mni_icbm152_nlin_asym_09c = "/path/to/atlas/folder2"
    )
  )

  # For T1 and postop CT coregistration only
  yael_preprocess(
    subject_code = "patient01",
    t1w_path = "/path/to/T1.nii or T1.nii.gz",
    ct_path = "/path/to/CT.nii or CT.nii.gz",

    # No normalization
    normalize_template = NULL,
    normalize_back = NA
  )

  # For both T1 and postop CT coregistration and T1 normalization
  yael_preprocess(
    subject_code = "patient01",
    t1w_path = "/path/to/T1.nii or T1.nii.gz",
    ct_path = "/path/to/CT.nii or CT.nii.gz",

    normalize_template = c(
      "mni_icbm152_nlin_asym_09a",
      "mni_icbm152_nlin_asym_09b",
      "mni_icbm152_nlin_asym_09c"
    ),

    normalize_policy = "auto",

    normalize_back = "mni_icbm152_nlin_asym_09b",

    atlases = list(
      mni_icbm152_nlin_asym_09b = "/path/to/atlas/folder1",
      mni_icbm152_nlin_asym_09c = "/path/to/atlas/folder2"
    )
  )

  ## End(Not run)

```

Description

Collapse power array with given analysis cubes

Usage

```
collapse_power(x, analysis_index_cubes)

## S3 method for class 'array'
collapse_power(x, analysis_index_cubes)

## S3 method for class 'FileArray'
collapse_power(x, analysis_index_cubes)
```

Arguments

x a [FileArray-class](#) array, must have 4 modes in the following sequence Frequency, Time, Trial, and Electrode

analysis_index_cubes a list of analysis indices for each mode

Value

a list of collapsed (mean) results

freq_trial_elec collapsed over time-points

freq_time_elec collapsed over trials

time_trial_elec collapsed over frequencies

freq_time collapsed over trials and electrodes

freq_elec collapsed over trials and time-points

freq_trial collapsed over time-points and electrodes

time_trial collapsed over frequencies and electrodes

time_elec collapsed over frequencies and trials

trial_elec collapsed over frequencies and time-points

freq power per frequency, averaged over other modes

time power per time-point, averaged over other modes

trial power per trial, averaged over other modes

Examples

```
if(!is_on_cran()) {

# Generate a 4-mode tensor array
x <- filearray::filearray_create(
  tempfile(), dimension = c(16, 100, 20, 5),
  partition_size = 1
)
```

```

x[] <- rnorm(160000)
dnames <- list(
  Frequency = 1:16,
  Time = seq(0, 1, length.out = 100),
  Trial = 1:20,
  Electrode = 1:5
)
dimnames(x) <- dnames

# Collapse array
results <- collapse_power(x, list(
  overall = list(),
  A = list(Trial = 1:5, Frequency = 1:6),
  B = list(Trial = 6:10, Time = 1:50)
))

# Plot power over frequency and time
groupB_result <- results$B

image(t(groupB_result$freq_time),
      x = dnames$Time[groupB_result$cube_index$Time],
      y = dnames$Frequency[groupB_result$cube_index$Frequency],
      xlab = "Time (s)",
      ylab = "Frequency (Hz)",
      xlim = range(dnames$Time))

x$delete(force = TRUE)

}

```

collapse2

*Collapse high-dimensional tensor array***Description**

Collapse high-dimensional tensor array

Usage

```

collapse2(x, keep, method = c("mean", "sum"), ...)

## S3 method for class 'FileArray'
collapse2(x, keep, method = c("mean", "sum"), ...)

## S3 method for class 'Tensor'
collapse2(x, keep, method = c("mean", "sum"), ...)

```

```
## S3 method for class 'array'
collapse2(x, keep, method = c("mean", "sum"), ...)
```

Arguments

x	R array, FileArray-class , or Tensor object
keep	integer vector, the margins to keep
method	character, calculates mean or sum of the array when collapsing
...	passed to other methods

Value

A collapsed array (or a vector or matrix), depending on keep

See Also

[collapse](#)

Examples

```
x <- array(1:16, rep(2, 4))

collapse2(x, c(3, 2))

# Alternative method, but slower when `x` is a large array
apply(x, c(3, 2), mean)

# filearray
y <- filearray::as_filearray(x)

collapse2(y, c(3, 2))

collapse2(y, c(3, 2), "sum")

# clean up
y$delete(force = TRUE)
```

compose_channel

Compose a "phantom" channel from existing electrodes

Description

In some cases, for example, deep-brain stimulation ('DBS'), it is often needed to analyze averaged electrode channels from segmented 'DBS' leads, or create bipolar contrast between electrode channels, or to generate non-equally weighted channel averages for 'Laplacian' reference. `compose_channel` allows users to generate a phantom channel that does not physically exist, but is treated as a normal electrode channel in 'RAVE'.

Usage

```
compose_channel(
  subject,
  number,
  from,
  weights = rep(1/length(from), length(from)),
  normalize = FALSE,
  force = FALSE,
  label = sprintf("Composed-%s", number),
  signal_type = c("auto", "LFP", "Spike", "EKG", "Auxiliary", "Unknown")
)
```

Arguments

subject	'RAVE' subject
number	new channel number, must be positive integer, cannot be existing electrode channel numbers
from	a vector of electrode channels that is used to compose this new channel, must be non-empty; see weights if these channels are not equally weighted.
weights	numerical weights used on each from channels; the length of weights must equals to the length of from; default is equally weighted for each channel (mean of from channels).
normalize	whether to normalize the weights such that the composed channel has the same variance as from channels; default is false
force	whether to overwrite existing composed channel if it exists; default is false. By specifying force=TRUE, users are risking breaking the data integrity since any analysis based on the composed channel is no longer reproducible. Also users cannot overwrite original channels under any circumstances.
label	the label for the composed channel; will be stored at 'electrodes.csv'
signal_type	signal type of the composed channel; default is 'auto' (same as the first from channel); other choices see SIGNAL_TYPES

Value

Nothing

Examples

```
library(raveio)

# Make sure demo subject exists in this example, just want to make
# sure the example does not error out
if(
  interactive() && "demo" %in% get_projects() &&
  "DemoSubject" %in% as_rave_project('demo')$subjects() &&
  local({
```

```

        subject <- as_rave_subject("demo/DemoSubject")
        !100 %in% subject$electrodes
    })
} {

# the actual example code:
# new channel 100 = 2 x channel 14 - (channe 15 + 16)
compose_channel(
  subject = "demo/DemoSubject",
  number = 100,
  from = c(14, 15, 16),
  weights = c(2, -1, -1),
  normalize = FALSE
)
}

```

convert_blackrock	<i>Convert 'BlackRock' 'NEV/NSx' files</i>
-------------------	--

Description

Convert 'BlackRock' 'NEV/NSx' files

Usage

```

convert_blackrock(
  file,
  block = NULL,
  subject = NULL,
  to = NULL,
  epoch = c("comment", "digital_inputs", "recording", "configuration", "log",
    "button_trigger", "tracking", "video_sync"),
  format = c("mat", "hdf5"),
  header_only = FALSE,
  ...
)

```

Arguments

file	path to any 'NEV/NSx' file
block	the block name, default is file name
subject	subject code to save the files; default is NULL

to	save to path, must be a directory; default is under the file path. If subject is provided, then the default is subject raw directory path
epoch	what type of events should be included in epoch file; default include comment, digital inputs, recording trigger, configuration change, log comment, button trigger, tracking, and video trigger.
format	output format, choices are 'mat' or 'hdf5'
header_only	whether just to generate channel and epoch table; default is false
...	ignored for enhanced backward compatibility

Value

The results will be stored in directory specified by to. Please read the output message carefully.

convert_electrode_table_to_bids
<i>Convert electrode table</i>

Description

Convert electrode table

Usage

```
convert_electrode_table_to_bids(  
  subject,  
  space = c("ScanRAS", "MNI305", "fsnative")  
)
```

Arguments

subject	'RAVE' subject
space	suggested coordinate space, notice this argument might not be supported when 'FreeSurfer' reconstruction is missing.

Value

A list of table in data frame and a list of meta information

dir_create2	<i>Force creating directory with checks</i>
-------------	---

Description

Force creating directory with checks

Usage

```
dir_create2(x, showWarnings = FALSE, recursive = TRUE, check = TRUE, ...)
```

Arguments

x	path to create
showWarnings, recursive, ...	passed to dir.create
check	whether to check the directory after creation

Value

Normalized path

Examples

```
path <- file.path(tempfile(), 'a', 'b', 'c')

# The following are equivalent
dir.create(path, showWarnings = FALSE, recursive = TRUE)

dir_create2(path)
```

ECoGTensor	<i>'iEEG/ECOG' Tensor class inherit from Tensor</i>
------------	---

Description

Four-mode tensor (array) especially designed for 'iEEG/ECOG' data. The Dimension names are: Trial, Frequency, Time, and Electrode.

Value

a data frame with the dimension names as index columns and value_name as value column
 an ECoGTensor instance

Super class

`raveio::Tensor` -> ECoGTensor

Methods**Public methods:**

- `ECoGTensor$flatten()`
- `ECoGTensor$new()`

Method `flatten()`: converts tensor (array) to a table (data frame)

Usage:

```
ECoGTensor$flatten(include_index = TRUE, value_name = "value")
```

Arguments:

`include_index` logical, whether to include dimension names

`value_name` character, column name of the value

Method `new()`: constructor

Usage:

```
ECoGTensor$new(
  data,
  dim,
  dimnames,
  varnames,
  hybrid = FALSE,
  swap_file = temp_tensor_file(),
  temporary = TRUE,
  multi_files = FALSE,
  use_index = TRUE,
  ...
)
```

Arguments:

`data` array or vector

`dim` dimension of data, must match with data

`dimnames` list of dimension names, equal length as `dim`

`varnames` names of `dimnames`, recommended names are: Trial, Frequency, Time, and Electrode

`hybrid` whether to enable hybrid mode to reduce RAM usage

`swap_file` if hybrid mode, where to store the data; default stores in `raveio_getopt('tensor_temp_path')`

`temporary` whether to clean up the space when exiting R session

`multi_files` logical, whether to use multiple files instead of one giant file to store data

`use_index` logical, when `multi_files` is true, whether use index dimension as partition number

`...` further passed to `Tensor` constructor

Author(s)

Zhengjia Wang

export_table

*Export data frame to different common formats***Description**

Stores and load data in various of data format. See 'Details' for limitations.

Usage

```
export_table(
  x,
  file,
  format = c("auto", "csv", "csv.zip", "h5", "fst", "json", "rds", "yaml"),
  ...
)

import_table(
  file,
  format = c("auto", "csv", "csv.zip", "h5", "fst", "json", "rds", "yaml"),
  ...
)
```

Arguments

x	data table to be saved to file
file	file to store the data
format	data storage format, default is 'auto' (infer from the file extension); other choices are 'csv', 'csv.zip', 'h5', 'fst', 'json', 'rds', 'yaml'
...	parameters passed to other functions

Details

The format 'rds', 'h5', 'fst', 'json', and 'yaml' try to preserve the first-level column attributes. Factors will be preserved in these formats. Such property does not exist in 'csv', 'csv.zip' formats.

Open-data formats are 'h5', 'csv', 'csv.zip', 'json', 'yaml'. These formats require the table elements to be native types (numeric, character, factor, etc.).

'rds', 'h5', and 'fst' can store large data sets. 'fst' is the best choice is performance and file size are the major concerns. 'rds' preserves all the properties of the table.

Value

The normalized path for export_table, and a [data.table](#) for import_table

Examples

```
x <- data.table::data.table(
  a = rnorm(10),
  b = letters[1:10],
  c = 1:10,
  d = factor(LETTERS[1:10])
)

f <- tempfile(fileext = ".csv.zip")

export_table(x = x, file = f)

y <- import_table(file = f)

str(x)
str(y)

# clean up
unlink(f)
```

find_path

Try to find path along the root directory

Description

Try to find path under root directory even if the original path is missing; see examples.

Usage

```
find_path(path, root_dir, all = FALSE)
```

Arguments

path	file path
root_dir	top directory of the search path
all	return all possible paths, default is false

Details

When file is missing, find_path concatenates the root directory and path combined to find the file. For example, if path is "a/b/c/d", the function first seek for existence of "a/b/c/d". If failed, then "b/c/d", and then "~/c/d" until reaching root directory. If all=TRUE, then all files/directories found along the search path will be returned

Value

The absolute path of file if exists, or NULL if missing/failed.

Examples

```

root <- tempdir()

# ----- Case 1: basic use case -----

# Create a path in root
dir_create2(file.path(root, 'a'))

# find path even it's missing. The search path will be
# root/ins/cd/a - missing
# root/cd/a      - missing
# root/a         - exists!
find_path('ins/cd/a', root)

# ----- Case 2: priority -----
# Create two paths in root
dir_create2(file.path(root, 'cc/a'))
dir_create2(file.path(root, 'a'))

# If two paths exist, return the first path found
# root/ins/cd/a - missing
# root/cd/a     - exists - returned
# root/a        - exists, but ignored
find_path('ins/cc/a', root)

# ----- Case 3: find all -----
# Create two paths in root
dir_create2(file.path(root, 'cc/a'))
dir_create2(file.path(root, 'a'))

# If two paths exist, return the first path found
# root/ins/cd/a - missing
# root/cd/a     - exists - returned
# root/a        - exists - returned
find_path('ins/cc/a', root, all = TRUE)

```

generate_reference	<i>Generate common average reference signal for 'RAVE' subjects</i>
--------------------	---

Description

To properly run this function, please install ravetools package.

Usage

```
generate_reference(subject, electrodes)
```

Arguments

subject	subject ID or RAVESubject instance
electrodes	electrodes to calculate the common average; these electrodes must run through 'Wavelet' first

Details

The goal of generating common average signals is to capture the common movement from all the channels and remove them out from electrode signals.

The common average signals will be stored at subject reference directories. Two exact same copies will be stored: one in 'HDF5' format such that the data can be read universally by other programming languages; one in [filearray](#) format that can be read in R with super fast speed.

Value

A reference instance returned by [new_reference](#) with signal type determined automatically.

get_projects	<i>Get all possible projects in 'RAVE' directory</i>
--------------	--

Description

Get all possible projects in 'RAVE' directory

Usage

```
get_projects(refresh = TRUE)
```

Arguments

refresh	whether to refresh the cache; default is true
---------	---

Value

characters of project names

get_val2	<i>Get value or return default if invalid</i>
----------	---

Description

Get value or return default if invalid

Usage

```
get_val2(x, key = NA, default = NULL, na = FALSE, min_len = 1L, ...)
```

Arguments

x	a list, or environment, or just any R object
key	the name to obtain from x. If NA, then return x. Default is NA
default	default value if
na, min_len, ...	passed to is_valid-ish

Value

values of the keys or default is invalid

Examples

```
x <- list(a=1, b = NA, c = character(0))

# ----- Basic usage -----

# no key, returns x if x is valid
get_val2(x)

get_val2(x, 'a', default = 'invalid')

# get 'b', NA is not filtered out
get_val2(x, 'b', default = 'invalid')

# get 'b', NA is considered invalid
get_val2(x, 'b', default = 'invalid', na = TRUE)

# get 'c', length 0 is allowed
get_val2(x, 'c', default = 'invalid', min_len = 0)

# length 0 is forbidden
get_val2(x, 'c', default = 'invalid', min_len = 1)
```

h5_names	<i>Returns all names contained in 'HDF5' file</i>
----------	---

Description

Returns all names contained in 'HDF5' file

Usage

```
h5_names(file)
```

Arguments

file	'HDF5' file path
------	------------------

Value

characters, data set names

h5_valid	<i>Check whether a 'HDF5' file can be opened for read/write</i>
----------	---

Description

Check whether a 'HDF5' file can be opened for read/write

Usage

```
h5_valid(file, mode = c("r", "w"), close_all = FALSE)
```

Arguments

file	path to file
mode	'r' for read access and 'w' for write access
close_all	whether to close all connections or just close current connection; default is false. Set this to TRUE if you want to close all other connections to the file

Value

logical whether the file can be opened.

Examples

```
x <- array(1:27, c(3,3,3))
f <- tempfile()

# No data written to the file, hence invalid
h5_valid(f, 'r')

save_h5(x, f, 'dset')
h5_valid(f, 'w') # TRUE
h5_valid(f, 'r') # TRUE
```

import_electrode_table

Import electrode table into subject meta folder

Description

Import electrode table into subject meta folder

Usage

```
import_electrode_table(path, subject, use_fs = NA, dry_run = FALSE, ...)
```

Arguments

path	path of table file, must be a 'csv' file
subject	'RAVE' subject ID or instance
use_fs	whether to use 'FreeSurfer' files to calculate other coordinates
dry_run	whether to dry-run the process; if true, then the table will be generated but not saved to subject's meta folder
...	passed to read.csv

Value

Nothing, the electrode information will be written directly to the subject's meta directory

ingest_regressor	<i>Ingest signals according to 'RAVE' repository epoch</i>
------------------	--

Description

Ingest signals according to 'RAVE' repository epoch

Usage

```
ingest_regressor(
  repository,
  source,
  source_type = c("channel", "file", "r_object"),
  filter = NULL,
  ...
)
```

Arguments

repository	'RAVE' repository, for example, generated from prepare_subject_power . The repository must contain epoch information.
source	data source, see source_type
source_type	type of the data source, choices are 'channel' source should be interpreted as subject channels (integer or series of integers) 'file' source should be interpreted as file path; the file format must be either 'Matlab' or 'HDF5', with data names to be the blocks in the "epoch" file. The sample rate of underlying signals must coincide with the sample rate presented in the repository 'r_object' source should be a list of R objects, where the names are blocks in the "epoch" file and the data are signal traces; The sample rate of underlying signals must coincide with the sample rate presented in the repository
filter	NULL or function; only used if source_type is 'channel'. For function, filter must have two exact arguments, with the first argument taking the signal data matrix (time point by channel) and second argument taking sample rates of length 2. The first sample rate is the original sampling frequency of the data; the second sample rate is sampling frequency derived from the repository.
...	passed to internal functions.

Value

A matrix of time by trial

install_subject	<i>Install a subject from the internet, a zip file or a directory</i>
-----------------	---

Description

Install a subject from the internet, a zip file or a directory

Usage

```
install_subject(
    path = ".",
    ask = interactive(),
    overwrite = FALSE,
    backup = TRUE,
    use_cache = TRUE,
    dry_run = FALSE,
    force_project = NA,
    force_subject = NA,
    ...
)
```

Arguments

path	path to subject archive, can be a path to directory, a zip file, or an internet address (must starts with 'http', or 'ftp')
ask	when overwrite is false, whether to ask the user if subject exists; default is true when running in interactive session; users will be prompt with choices; if ask=FALSE and overwrite=FALSE, then the process will end with a warning if the subject exists.
overwrite	whether to overwrite existing subject, see argument ask and backup
backup	whether to back-up the subject when overwriting the data; default is true, which will rename the old subject folders instead of removing; set to true to remove existing subject.
use_cache	whether to use cached extraction directory; default is true. Set it to FALSE if you want a clean installation.
dry_run	whether to dry-run the process instead of actually installing; this rehearsal can help you see the progress and prevent you from losing data
force_project, force_subject	force set the project or subject; will raise a warning as this might mess up some pipelines
...	passed to download.file

Examples

```
# Please run 2nd example of function archive_subject

## Not run:

install_subject(path)

## End(Not run)
```

is.blank	<i>Check If Input Has Blank String</i>
----------	--

Description

Check If Input Has Blank String

Usage

```
is.blank(x)
```

Arguments

x input data: a vector or an array

Value

```
x == ""
```

is.zerolenth	<i>Check If Input Has Zero Length</i>
--------------	---------------------------------------

Description

Check If Input Has Zero Length

Usage

```
is.zerolenth(x)
```

Arguments

x input data: a vector, list, or array

Value

whether x has zero length

is_on_cran	<i>Check if current session is on 'CRAN'</i>
------------	--

Description

Use this function only for examples and test. The goal is to comply with the 'CRAN' policy. Do not use it in normal functions to cheat. Violating 'CRAN' policy will introduce instability to your code. Make sure reading Section 'Details' before using this function.

Usage

```
is_on_cran(if_interactive = FALSE, verbose = FALSE)
```

Arguments

`if_interactive` whether interactive session will be considered as on 'CRAN'; default is FALSE
`verbose` whether to print out reason of return; default is no

Details

According to 'CRAN' policy, package examples and test functions may only use maximum 2 'CPU' cores. Examples running too long should be suppressed. Normally package developers will use `interactive()` to avoid running examples or parallel code on 'CRAN'. However, when checked locally, these examples will be skipped too. Coding bug in those examples will not be reported.

The objective is to allow 'RAVE' package developers to write and test examples locally or on integrated development environment (such as 'Github'), while suppressing them on 'CRAN'. In such way, bugs in the examples will be revealed and fixed promptly.

Do not use this function inside of the package functions to cheat or slip illegal code under the eyes of 'CRAN' folks. This will increase their work load and introduce instability to your code. Also this function is subject to deletion in the future.

Value

A logical whether current environment should be considered as on 'CRAN'.

is_valid_ish	<i>Check if data is close to "valid"</i>
--------------	--

Description

Check if data is close to "valid"

Usage

```
is_valid_ish(
  x,
  min_len = 1,
  max_len = Inf,
  mode = NA,
  na = TRUE,
  blank = FALSE,
  all = FALSE
)
```

Arguments

x	data to check
min_len, max_len	minimal and maximum length
mode	which storage mode (see mode) should x be considered valid. Default is NA: disabled.
na	whether NA values considered invalid?
blank	whether blank string considered invalid?
all	if na or blank is true, whether all element of x being invalid will result in failure?

Value

logicals whether input x is valid

Examples

```
# length checks
is_valid_ish(NULL) # FALSE
is_valid_ish(integer(0)) # FALSE
is_valid_ish(integer(0), min_len = 0) # TRUE
is_valid_ish(1:10, max_len = 9) # FALSE

# mode check
is_valid_ish(1:10) # TRUE
is_valid_ish(1:10, mode = 'numeric') # TRUE
is_valid_ish(1:10, mode = 'character') # FALSE

# NA or blank checks
is_valid_ish(NA) # FALSE
is_valid_ish(c(1,2,NA), all = FALSE) # FALSE
is_valid_ish(c(1,2,NA), all = TRUE) # TRUE as not all elements are NA

is_valid_ish(c('1',''), all = FALSE) # TRUE
is_valid_ish(1:3, all = FALSE) # TRUE as 1:3 are not characters
```

join_tensors*Join Multiple Tensors into One Tensor*

Description

Join Multiple Tensors into One Tensor

Usage

```
join_tensors(tensors, temporary = TRUE)
```

Arguments

tensors	list of Tensor instances
temporary	whether to garbage collect space when exiting R session

Details

Merges multiple tensors. Each tensor must share the same dimension with the last one dimension as 1, for example, 100x100x1. Join 3 tensors like this will result in a 100x100x3 tensor. This function is handy when each sub-tensors are generated separately. However, it does no validation test. Use with cautions.

Value

A new [Tensor](#) instance with the last dimension

Author(s)

Zhengjia Wang

Examples

```
tensor1 <- Tensor$new(data = 1:9, c(3,3,1), dimnames = list(
  A = 1:3, B = 1:3, C = 1
), varnames = c('A', 'B', 'C'))
tensor2 <- Tensor$new(data = 10:18, c(3,3,1), dimnames = list(
  A = 1:3, B = 1:3, C = 2
), varnames = c('A', 'B', 'C'))
merged <- join_tensors(list(tensor1, tensor2))
merged$get_data()
```

lapply_async

Run [lapply](#) in parallel

Description

Uses [lapply_async2](#), but allows better parallel scheduling via [with_future_parallel](#). On 'Unix', the function will fork processes. On 'Windows', the function uses strategies specified by `on_failure`

Usage

```
lapply_async(
  x,
  FUN,
  FUN.args = list(),
  callback = NULL,
  ncores = NULL,
  on_failure = "multisession",
  ...
)
```

Arguments

<code>x</code>	iterative elements
<code>FUN</code>	function to apply to each element of <code>x</code>
<code>FUN.args</code>	named list that will be passed to <code>FUN</code> as arguments
<code>callback</code>	callback function or <code>NULL</code> . When passed as function, the function takes one argument (elements of <code>x</code>) as input, and it suppose to return one string character.
<code>ncores</code>	number of cores to use, constraint by the <code>max_worker</code> option (see raveio_getopt); default is the maximum number of workers available
<code>on_failure</code>	alternative strategy if fork process is disallowed (set by users or on 'Windows')
<code>...</code>	passed to lapply_async2

Examples

```
if(!is_on_cran()) {
  library(raveio)

  # ---- Basic example -----
  lapply_async(1:16, function(x) {
    # function that takes long to fun
    Sys.sleep(1)
    x
  })

  # With callback
```

```

lapply_async(1:16, function(x){
  Sys.sleep(1)
  x + 1
}, callback = function(x) {
  sprintf("Calculating|%", x)
})

# With ncores
pids <- lapply_async(1:16, function(x){
  Sys.sleep(0.5)
  Sys.getpid()
}, ncores = 2)

# Unique number of PIDs (cores)
unique(unlist(pids))

# ---- With scheduler -----
# Scheduler pre-initialize parallel workers and temporary
# switches parallel context. The workers ramp-up
# time can be saved by reusing the workers.
#
with_future_parallel({

  # lapply_async block 1
  pids <- lapply_async(1:16, function(x){
    Sys.sleep(1)
    Sys.getpid()
  }, callback = function(x) {
    sprintf("lapply_async without ncores|%", x)
  })
  print(unique(unlist(pids)))

  # lapply_async block 2
  pids <- lapply_async(1:16, function(x){
    Sys.sleep(1)
    Sys.getpid()
  }, callback = function(x) {
    sprintf("lapply_async with ncores|%", x)
  }, ncores = 4)
  print(unique(unlist(pids)))

})

}

```

Description

Please use a safer [new_electrode](#) function to create instances. This documentation is to describe the member methods of the electrode class LFP_electrode

Value

if the reference number is NULL or 'noref', then returns 0, otherwise returns a [FileArray-class](#)

If simplify is enabled, and only one block is loaded, then the result will be a vector (type="voltage") or a matrix (others), otherwise the result will be a named list where the names are the blocks.

Super class

[raveio::RAVEAbstractElectrode](#) -> LFP_electrode

Active bindings

h5_fname 'HDF5' file name

valid whether current electrode is valid: subject exists and contains current electrode or reference;
subject electrode type matches with current electrode type

raw_sample_rate voltage sample rate

power_sample_rate power/phase sample rate

preprocess_info preprocess information

power_file path to power 'HDF5' file

phase_file path to phase 'HDF5' file

voltage_file path to voltage 'HDF5' file

Methods**Public methods:**

- [LFP_electrode\\$print\(\)](#)
- [LFP_electrode\\$set_reference\(\)](#)
- [LFP_electrode\\$new\(\)](#)
- [LFP_electrode\\$.load_noref_wavelet\(\)](#)
- [LFP_electrode\\$.load_noref_voltage\(\)](#)
- [LFP_electrode\\$.load_wavelet\(\)](#)
- [LFP_electrode\\$.load_voltage\(\)](#)
- [LFP_electrode\\$.load_raw_voltage\(\)](#)
- [LFP_electrode\\$load_data\(\)](#)
- [LFP_electrode\\$load_blocks\(\)](#)
- [LFP_electrode\\$clear_cache\(\)](#)
- [LFP_electrode\\$clear_memory\(\)](#)
- [LFP_electrode\\$clone\(\)](#)

Method [print\(\)](#): print electrode summary

Usage:

```
LFP_electrode$print()
```

Method `set_reference()`: set reference for current electrode

Usage:

```
LFP_electrode$set_reference(reference)
```

Arguments:

reference either NULL or LFP_electrode instance

Method `new()`: constructor

Usage:

```
LFP_electrode$new(subject, number, quiet = FALSE)
```

Arguments:

subject, number, quiet see constructor in [RAVEAbstarctElectrode](#)

Method `.load_noref_wavelet()`: load non-referenced wavelet coefficients (internally used)

Usage:

```
LFP_electrode$.load_noref_wavelet(reload = FALSE)
```

Arguments:

reload whether to reload cache

Method `.load_noref_voltage()`: load non-referenced voltage (internally used)

Usage:

```
LFP_electrode$.load_noref_voltage(reload = FALSE)
```

Arguments:

reload whether to reload cache

srate voltage signal sample rate

Method `.load_wavelet()`: load referenced wavelet coefficients (internally used)

Usage:

```
LFP_electrode$.load_wavelet(  
  type = c("power", "phase", "wavelet-coefficient"),  
  reload = FALSE  
)
```

Arguments:

type type of data to load

reload whether to reload cache

Method `.load_voltage()`: load referenced voltage (internally used)

Usage:

```
LFP_electrode$.load_voltage(reload = FALSE)
```

Arguments:

reload whether to reload cache

Method `.load_raw_voltage()`: load raw voltage (no process)

Usage:

```
LFP_electrode$.load_raw_voltage(reload = FALSE)
```

Arguments:

`reload` whether to reload cache

Method `load_data()`: method to load electrode data

Usage:

```
LFP_electrode$load_data(
  type = c("power", "phase", "voltage", "wavelet-coefficient", "raw-voltage")
)
```

Arguments:

`type` data type such as "power", "phase", "voltage", "wavelet-coefficient", and "raw-voltage".

For "power", "phase", and "wavelet-coefficient", 'Wavelet' transforms are required.

For "voltage", 'Notch' filters must be applied. All these types except for "raw-voltage" will be referenced. For "raw-voltage", no reference will be performed since the data will be the "raw" signal (no processing).

Method `load_blocks()`: load electrode block-wise data (with no reference), useful when epoch is absent

Usage:

```
LFP_electrode$load_blocks(
  blocks,
  type = c("power", "phase", "voltage", "wavelet-coefficient", "raw-voltage"),
  simplify = TRUE
)
```

Arguments:

`blocks` session blocks

`type` data type such as "power", "phase", "voltage", "raw-voltage" (with no filters applied, as-is from imported), "wavelet-coefficient". Note that if type is "raw-voltage", then the data only needs to be imported; for "voltage" data, 'Notch' filters must be applied; for all other types, 'Wavelet' transforms are required.

`simplify` whether to simplify the result

Method `clear_cache()`: method to clear cache on hard drive

Usage:

```
LFP_electrode$clear_cache(...)
```

Arguments:

`...` ignored

Method `clear_memory()`: method to clear memory

Usage:

```
LFP_electrode$clear_memory(...)
```

Arguments:

... ignored

Method clone(): The objects of this class are cloneable with this method.

Usage:

```
LFP_electrode$clone(deep = FALSE)
```

Arguments:

deep Whether to make a deep clone.

Examples

```
# Download subject demo/DemoSubject

subject <- as_rave_subject("demo/DemoSubject", strict = FALSE)

if(dir.exists(subject$path)) {

# Electrode 14 in demo/DemoSubject
e <- new_electrode(subject = subject, number = 14, signal_type = "LFP")

# Load CAR reference "ref_13-16,24"
ref <- new_reference(subject = subject, number = "ref_13-16,24",
                     signal_type = "LFP")
e$set_reference(ref)

# Set epoch
e$set_epoch(epoch = 'auditory_onset')

# Set loading window
e$trial_intervals <- list(c(-1, 2))

# Preview
print(e)

# Now epoch power
power <- e$load_data("power")
names(dimnames(power))

# Subset power
subset(power, Time ~ Time < 0, Electrode ~ Electrode == 14)

# Draw baseline
tempfile <- tempfile()
bl <- power_baseline(power, baseline_windows = c(-1, 0),
                    method = "decibel", filebase = tempfile)
collapsed_power <- collapse2(bl, keep = c(2,1))
# Visualize
dname <- dimnames(bl)
image(collapsed_power, x = dname$Time, y = dname$Frequency,
      xlab = "Time (s)", ylab = "Frequency (Hz)",
      main = "Mean power over trial (Baseline: -1~0 seconds)",
      sub = glue('Electrode {e$number} (Reference: {ref$number}')))
```

```

abline(v = 0, lty = 2, col = 'blue')
text(x = 0, y = 20, "Audio onset", col = "blue", cex = 0.6)

# clear cache on hard disk
e$clear_cache()
ref$clear_cache()

}

```

LFP_reference

Definitions of reference with 'LFP' signal type

Description

Please use a safer [new_reference](#) function to create instances. This documentation is to describe the member methods of the electrode class LFP_reference

Value

if the reference number is NULL or 'noref', then returns 0, otherwise returns a [FileArray-class](#)

If simplify is enabled, and only one block is loaded, then the result will be a vector (type="voltage") or a matrix (others), otherwise the result will be a named list where the names are the blocks.

Super class

[raveio::RAVEAbstractElectrode](#) -> LFP_reference

Active bindings

exists whether electrode exists in subject

h5_fname 'HDF5' file name

valid whether current electrode is valid: subject exists and contains current electrode or reference;
subject electrode type matches with current electrode type

raw_sample_rate voltage sample rate

power_sample_rate power/phase sample rate

preprocess_info preprocess information

power_file path to power 'HDF5' file

phase_file path to phase 'HDF5' file

voltage_file path to voltage 'HDF5' file

Methods

Public methods:

- `LFP_reference$print()`
- `LFP_reference$set_reference()`
- `LFP_reference$new()`
- `LFP_reference$.load_noref_wavelet()`
- `LFP_reference$.load_noref_voltage()`
- `LFP_reference$.load_wavelet()`
- `LFP_reference$.load_voltage()`
- `LFP_reference$load_data()`
- `LFP_reference$load_blocks()`
- `LFP_reference$clear_cache()`
- `LFP_reference$clear_memory()`
- `LFP_reference$clone()`

Method `print()`: print reference summary

Usage:

`LFP_reference$print()`

Method `set_reference()`: set reference for current electrode

Usage:

`LFP_reference$set_reference(reference)`

Arguments:

reference either NULL or LFP_electrode instance

Method `new()`: constructor

Usage:

`LFP_reference$new(subject, number, quiet = FALSE)`

Arguments:

subject, number, quiet see constructor in [RAVEAbstarctElectrode](#)

Method `.load_noref_wavelet()`: load non-referenced wavelet coefficients (internally used)

Usage:

`LFP_reference$.load_noref_wavelet(reload = FALSE)`

Arguments:

reload whether to reload cache

Method `.load_noref_voltage()`: load non-referenced voltage (internally used)

Usage:

`LFP_reference$.load_noref_voltage(reload = FALSE)`

Arguments:

reload whether to reload cache

srate voltage signal sample rate

Method `.load_wavelet()`: load referenced wavelet coefficients (internally used)

Usage:

```
LFP_reference$.load_wavelet(
  type = c("power", "phase", "wavelet-coefficient"),
  reload = FALSE
)
```

Arguments:

type type of data to load
reload whether to reload cache

Method `.load_voltage()`: load referenced voltage (internally used)

Usage:

```
LFP_reference$.load_voltage(reload = FALSE)
```

Arguments:

reload whether to reload cache

Method `load_data()`: method to load electrode data

Usage:

```
LFP_reference$load_data(
  type = c("power", "phase", "voltage", "wavelet-coefficient")
)
```

Arguments:

type data type such as "power", "phase", "voltage", "wavelet-coefficient".

Method `load_blocks()`: load electrode block-wise data (with reference), useful when epoch is absent

Usage:

```
LFP_reference$load_blocks(
  blocks,
  type = c("power", "phase", "voltage", "wavelet-coefficient"),
  simplify = TRUE
)
```

Arguments:

blocks session blocks
type data type such as "power", "phase", "voltage", "wavelet-coefficient". Note that if type is voltage, then 'Notch' filters must be applied; otherwise 'Wavelet' transforms are required.
simplify whether to simplify the result

Method `clear_cache()`: method to clear cache on hard drive

Usage:

```
LFP_reference$clear_cache(...)
```

Arguments:

... ignored

Method `clear_memory()`: method to clear memory

Usage:

`LFP_reference$clear_memory(...)`

Arguments:

... ignored

Method `clone()`: The objects of this class are cloneable with this method.

Usage:

`LFP_reference$clone(deep = FALSE)`

Arguments:

`deep` Whether to make a deep clone.

Examples

```
## Not run:

# Download subject demo/DemoSubject

subject <- as_rave_subject("demo/DemoSubject")

# Electrode 14 as reference electrode (Bipolar referencing)
e <- new_reference(subject = subject, number = "ref_14",
  signal_type = "LFP")

# Reference "ref_13-16,24" (CAR or white-matter reference)
ref <- new_reference(subject = subject, number = "ref_13-16,24",
  signal_type = "LFP")
ref

# Set epoch
e$set_epoch(epoch = 'auditory_onset')

# Set loading window
e$trial_intervals <- list(c(-1, 2))

# Preview
print(e)

# Now epoch power
power <- e$load_data("power")
names(dimnames(power))

# Subset power
subset(power, Time ~ Time < 0, Electrode ~ Electrode == 14)

# clear cache on hard disk
```

```
e$clear_cache()
```

```
## End(Not run)
```

load_bids_ieeg_header *Read in description files from 'BIDS-iEEG' format*

Description

Analyze file structures and import all json and tsv files. File specification can be found at <https://bids-specification.readthedocs.io/en/stable/>, chapter "Modality specific files", section "Intracranial Electroencephalography" ([doi:10.1038/s4159701901057](https://doi.org/10.1038/s4159701901057)). Please note that this function has very limited support on BIDS format.

Usage

```
load_bids_ieeg_header(bids_root, project_name, subject_code, folder = "ieeg")
```

Arguments

bids_root	'BIDS' root directory
project_name	project folder name
subject_code	subject code, do not include "sub-" prefix
folder	folder name corresponding to 'iEEG' data. It's possible to analyze other folders. However, by default, the function is designed for 'ieeg' folder.

Value

A list containing the information below:

subject_code	character, removed leading "sub-"
project_name	character, project name
has_session	whether session/block names are indicated by the file structure
session_names	session/block names indicated by file structure. If missing, then session name will be "default"
paths	a list containing path information
stimuli_path	stimuli path, not used for now
sessions	A named list containing meta information for each session/block. The names of the list is task name, and the items corresponding to the task contains events and channel information. Miscellaneous files are stored in "others" variable.

Examples

```
# Download https://github.com/bids-standard/bids-examples/
# extract to directory ~/rave_data/bids_dir/

bids_root <- '~/rave_data/bids_dir/'
project_name <- 'ieeg_visual'

if(dir.exists(bids_root) &&
  dir.exists(file.path(bids_root, project_name, 'sub-01'))){

  header <- load_bids_ieeg_header(bids_root, project_name, '01')

  print(header)

  # sessions
  names(header$sessions)

  # electrodes
  head(header$sessions$`01`$spaces$unknown_space$table)

  # visual task channel settings
  head(header$sessions$`01`$tasks$`01-visual-01`$channels)

  # event table
  head(header$sessions$`01`$tasks$`01-visual-01`$channels)
}
```

load_fst_or_h5

*Function try to load 'fst' arrays, if not found, read 'HDF5' arrays***Description**

Function try to load 'fst' arrays, if not found, read 'HDF5' arrays

Usage

```
load_fst_or_h5(
  fst_path,
  h5_path,
  h5_name,
  fst_need_transpose = FALSE,
  fst_need_drop = FALSE,
  ram = FALSE
)
```

Arguments

fst_path	'fst' file cache path
h5_path	alternative 'HDF5' file path
h5_name	'HDF5' data name
fst_need_transpose	does 'fst' data need transpose?
fst_need_drop	drop dimensions
ram	whether to load to memory directly or perform lazy loading

Details

RAVE stores data with redundancy. One electrode data is usually saved with two copies in different formats: 'HDF5' and 'fst', where 'HDF5' is cross-platform and supported by multiple languages such as Matlab, Python, etc, while 'fst' format is supported by R only, with super high read/write speed. `load_fst_or_h5` checks whether the presence of 'fst' file, if failed, then it reads data from persistent 'HDF5' file.

Value

If 'fst' cache file exists, returns [LazyFST](#) object, otherwise returns [LazyH5](#) instance

load_h5	<i>Lazy Load 'HDF5' File via io_read_h5</i>
---------	---

Description

Wrapper for class [LazyH5](#), which load data with "lazy" mode - only read part of dataset when needed.

Usage

```
load_h5(file, name, read_only = TRUE, ram = FALSE, quiet = FALSE)
```

Arguments

file	'HDF5' file
name	group/data_name path to dataset (H5D data)
read_only	only used if ram=FALSE, whether the returned LazyH5 instance should be read only
ram	load data to memory immediately, default is false
quiet	whether to suppress messages

Value

If ram is true, then return data as arrays, otherwise return a [LazyH5](#) instance.

See Also[save_h5](#)**Examples**

```

file <- tempfile()
x <- array(1:120, dim = c(4,5,6))

# save x to file with name /group/dataset/1
save_h5(x, file, '/group/dataset/1', quiet = TRUE)

# read data
y <- load_h5(file, '/group/dataset/1', ram = TRUE)
class(y) # array

z <- load_h5(file, '/group/dataset/1', ram = FALSE)
class(z) # LazyH5

dim(z)

```

load_meta2

*Load 'RAVE' subject meta data***Description**

Load 'RAVE' subject meta data

Usage

```
load_meta2(meta_type, project_name, subject_code, subject_id, meta_name)
```

Arguments

meta_type	electrodes, epochs, time_points, frequencies, references ...
project_name	project name
subject_code	subject code
subject_id	"project_name/subject_code"
meta_name	only used if meta_type is epochs or references

Value

A data frame of the specified meta type or NULL is no meta data is found.

load_yaml	A port to read_yaml
-----------	-------------------------------------

Description

For more examples, see [save_yaml](#).

Usage

```
load_yaml(file, ..., map = NULL)
```

Arguments

file, ...	passed to read_yaml
map	fastmap2 instance or NULL

Value

A [fastmap2](#). If map is provided then return map, otherwise return newly created one

See Also

[fastmap2](#), [save_yaml](#), [read_yaml](#), [write_yaml](#)

mgf_to_nii	Convert 'FreeSurfer' 'mgf' to 'Nifti'
------------	---------------------------------------

Description

Convert 'FreeSurfer' 'mgf' to 'Nifti'

Usage

```
mgf_to_nii(from, to)
```

Arguments

from	path to 'FreeSurfer' 'mgf' or 'mgz' file
to	path to 'Nifti' file, must ends with 'nii' or 'nii.gz'

Value

Nothing; the file will be created to path specified by to

new_electrode	Create new electrode channel instance or a reference signal instance
---------------	--

Description

Create new electrode channel instance or a reference signal instance

Usage

```
new_electrode(subject, number, signal_type, ...)
```

```
new_reference(subject, number, signal_type, ...)
```

Arguments

subject	characters, or a RAVESubject instance
number	integer in new_electrode, or characters in new_reference; see 'Details' and 'Examples'
signal_type	signal type of the electrode or reference; can be automatically inferred, but it is highly recommended to specify a value; see SIGNAL_TYPES
...	other parameters passed to class constructors, respectively

Details

In new_electrode, number should be a positive valid integer indicating the electrode number. In new_reference, number can be one of the followings:

'noref', or NULL no reference is needed

'ref_X' where 'X' is a single number, then the reference is another existing electrode; this could occur in bipolar-reference cases

'ref_XXX' 'XXX' is a combination of multiple electrodes that can be parsed by [parse_svec](#). This could occur in common average reference, or white matter reference. One example is 'ref_13-16,24', meaning the reference signal is an average of electrode 13, 14, 15, 16, and 24.

Value

Electrode or reference instances that inherit [RAVEAbstarctElectrode](#) class

Examples

```
## Not run:

# Download subject demo/DemoSubject (~500 MB)

# Electrode 14 in demo/DemoSubject
subject <- as_rave_subject("demo/DemoSubject")
```

```

e <- new_electrode(subject = subject, number = 14, signal_type = "LFP")

# Load CAR reference "ref_13-16,24"
ref <- new_reference(subject = subject, number = "ref_13-16,24",
                     signal_type = "LFP")
e$set_reference(ref)

# Set epoch
e$set_epoch(epoch = 'auditory_onset')

# Set loading window
e$trial_intervals <- list(c(-1, 2))

# Preview
print(e)

# Now epoch power
power <- e$load_data("power")
names(dimnames(power))

# Subset power
subset(power, Time ~ Time < 0, Electrode ~ Electrode == 14)

# Draw baseline
tempfile <- tempfile()
bl <- power_baseline(power, baseline_windows = c(-1, 0),
                     method = "decibel", filebase = tempfile)
collapsed_power <- collapse2(bl, keep = c(2,1))
# Visualize
dname <- dimnames(bl)
image(collapsed_power, x = dname$Time, y = dname$Frequency,
      xlab = "Time (s)", ylab = "Frequency (Hz)",
      main = "Mean power over trial (Baseline: -1~0 seconds)",
      sub = glue('Electrode {e$number} (Reference: {ref$number})'))
abline(v = 0, lty = 2, col = 'blue')
text(x = 0, y = 20, "Audio onset", col = "blue", cex = 0.6)

# clear cache on hard disk
e$clear_cache()
ref$clear_cache()

## End(Not run)

```

niftyreg_coreg

Register 'CT' to 'MR' images via 'NiftyReg'

Description

Supports 'Rigid', 'affine', or 'non-linear' transformation

Usage

```

niftyreg_coreg(
  ct_path,
  mri_path,
  coreg_path = NULL,
  reg_type = c("rigid", "affine", "nonlinear"),
  interp = c("trilinear", "cubic", "nearest"),
  verbose = TRUE,
  ...
)

cmd_run_niftyreg_coreg(
  subject,
  ct_path,
  mri_path,
  reg_type = c("rigid", "affine", "nonlinear"),
  interp = c("trilinear", "cubic", "nearest"),
  verbose = TRUE,
  dry_run = FALSE,
  ...
)

```

Arguments

<code>ct_path, mri_path</code>	absolute paths to 'CT' and 'MR' image files
<code>coreg_path</code>	registration path, where to save results; default is the parent folder of <code>ct_path</code>
<code>reg_type</code>	registration type, choices are 'rigid', 'affine', or 'nonlinear'
<code>interp</code>	how to interpolate when sampling volumes, choices are 'trilinear', 'cubic', or 'nearest'
<code>verbose</code>	whether to verbose command; default is true
<code>...</code>	other arguments passed to register_volume
<code>subject</code>	'RAVE' subject
<code>dry_run</code>	whether to dry-run the script and to print out the command instead of executing the code; default is false

Value

Nothing is returned from the function. However, several files will be generated at the 'CT' path:

- 'ct_in_t1.nii' aligned 'CT' image; the image is also re-sampled into 'MRI' space
- 'CT_IJK_to_MR_RAS.txt' transform matrix from volume 'IJK' space in the original 'CT' to the 'RAS' anatomical coordinate in 'MR' scanner
- 'CT_RAS_to_MR_RAS.txt' transform matrix from scanner 'RAS' space in the original 'CT' to 'RAS' in 'MR' scanner space

power_baseline	<i>Calculate power baseline</i>
----------------	---------------------------------

Description

Calculate power baseline

Usage

```
power_baseline(
  x,
  baseline_windows,
  method = c("percentage", "sqrt_percentage", "decibel", "zscore", "sqrt_zscore"),
  units = c("Trial", "Frequency", "Electrode"),
  ...
)

## S3 method for class 'rave_prepare_power'
power_baseline(
  x,
  baseline_windows,
  method = c("percentage", "sqrt_percentage", "decibel", "zscore", "sqrt_zscore"),
  units = c("Frequency", "Trial", "Electrode"),
  electrodes,
  ...
)

## S3 method for class 'FileArray'
power_baseline(
  x,
  baseline_windows,
  method = c("percentage", "sqrt_percentage", "decibel", "zscore", "sqrt_zscore"),
  units = c("Frequency", "Trial", "Electrode"),
  filebase = NULL,
  ...
)

## S3 method for class 'array'
power_baseline(
  x,
  baseline_windows,
  method = c("percentage", "sqrt_percentage", "decibel", "zscore", "sqrt_zscore"),
  units = c("Trial", "Frequency", "Electrode"),
  ...
)

## S3 method for class 'ECoGTensor'
```

```

power_baseline(
  x,
  baseline_windows,
  method = c("percentage", "sqrt_percentage", "decibel", "zscore", "sqrt_zscore"),
  units = c("Trial", "Frequency", "Electrode"),
  filebase = NULL,
  hybrid = TRUE,
  ...
)

```

Arguments

<code>x</code>	R array, filearray , ECoGTensor , or 'rave_prepare_power' object created by prepare_subject_power .
<code>baseline_windows</code>	list of baseline window (intervals)
<code>method</code>	baseline method; choices are 'percentage', 'sqrt_percentage', 'decibel', 'zscore', 'sqrt_zscore'; see 'Details' in baseline_array
<code>units</code>	the unit of the baseline; see 'Details'
<code>...</code>	passed to other methods
<code>electrodes</code>	the electrodes to be included in baseline calculation; for power repository object produced by prepare_subject_power only; default is all available electrodes in each of <code>signal_types</code>
<code>filebase</code>	where to store the output; default is NULL and is automatically determined
<code>hybrid</code>	whether the array will be

Details

The arrays must be four-mode tensor and must have valid named [dimnames](#). The dimension names must be 'Trial', 'Frequency', 'Time', 'Electrode', case sensitive.

The `baseline_windows` determines the baseline windows that are used to calculate time-points of baseline to be included. This can be one or more intervals and must pass the validation function [validate_time_window](#).

The `units` determines the unit of the baseline. It can be one or more of 'Trial', 'Frequency', 'Electrode'. The default value is all of them, i.e., baseline for each combination of trial, frequency, and electrode. To share the baseline across trials, please remove 'Trial' from units. To calculate baseline that should be shared across electrodes (e.g. in some mini-electrodes), remove 'Electrode' from the units.

Value

Usually the same type as the input: for arrays, [filearray](#), or [ECoGTensor](#), the outputs are also the same type with the same dimensions; for 'rave_prepare_power' repositories, the results will be stored in its 'baselined' element; see 'Examples'.

Examples

```
## Not run:
# The following code need to download additional demo data
# Please see https://rave.wiki/ for more details

library(raveio)
repo <- prepare_subject_power(
  subject = "demo/DemoSubject",
  time_windows = c(-1, 3),
  electrodes = c(14, 15))

##### Direct baseline on the repository
power_baseline(x = repo, method = "decibel",
               baseline_windows = list(c(-1, 0), c(2, 3)))
power_mean <- repo$power$baselined$collapse(
  keep = c(2,1), method = "mean")
image(power_mean, x = repo$time_points, y = repo$frequency,
       xlab = "Time (s)", ylab = "Frequency (Hz)",
       main = "Mean power over trial (Baseline: -1~0 & 2~3)")
abline(v = 0, lty = 2, col = 'blue')
text(x = 0, y = 20, "Aud-Onset", col = "blue", cex = 0.6)

##### Alternatively, baseline on electrode instances
baselined <- lapply(repo$power$data_list, function(inst) {
  re <- power_baseline(inst, method = "decibel",
                       baseline_windows = list(c(-1, 0), c(2, 3)))
  collapse2(re, keep = c(2,1), method = "mean")
})
power_mean2 <- (baselined[[1]] + baselined[[2]]) / 2

# Same with precision difference
max(abs(power_mean2 - power_mean)) < 1e-6

## End(Not run)
```

prepare_subject_bare0 *Prepare 'RAVE' single-subject data*

Description

Prepare 'RAVE' single-subject data

Usage

```
prepare_subject_bare0(
  subject,
```

```
    electrodes,  
    reference_name,  
    ...,  
    quiet = TRUE,  
    repository_id = NULL  
)
```

```
prepare_subject_bare(  
  subject,  
  electrodes,  
  reference_name,  
  ...,  
  repository_id = NULL  
)
```

```
prepare_subject_with_epoch(  
  subject,  
  electrodes,  
  reference_name,  
  epoch_name,  
  time_windows,  
  stitch_events = NULL,  
  env = parent.frame(),  
  ...  
)
```

```
prepare_subject_with_blocks(  
  subject,  
  electrodes,  
  reference_name,  
  blocks,  
  raw = FALSE,  
  signal_type = "LFP",  
  time_frequency = (!raw && signal_type == "LFP"),  
  quiet = raw,  
  env = parent.frame(),  
  repository_id = NULL,  
  ...  
)
```

```
prepare_subject_phase(  
  subject,  
  electrodes,  
  reference_name,  
  epoch_name,  
  time_windows,  
  stitch_events = NULL,  
  signal_type = c("LFP"),
```

```
env = parent.frame(),
verbose = TRUE,
...
)

prepare_subject_power(
  subject,
  electrodes,
  reference_name,
  epoch_name,
  time_windows,
  stitch_events = NULL,
  signal_type = c("LFP"),
  env = parent.frame(),
  verbose = TRUE,
  ...
)

prepare_subject_wavelet(
  subject,
  electrodes,
  reference_name,
  epoch_name,
  time_windows,
  stitch_events = NULL,
  signal_type = c("LFP"),
  env = parent.frame(),
  verbose = TRUE,
  ...
)

prepare_subject_raw_voltage_with_epoch(
  subject,
  electrodes,
  epoch_name,
  time_windows,
  stitch_events = NULL,
  ...,
  quiet = TRUE,
  repository_id = NULL
)

prepare_subject_voltage_with_epoch(
  subject,
  electrodes,
  epoch_name,
  time_windows,
  reference_name,
```

```

    stitch_events = NULL,
    ...,
    quiet = TRUE,
    repository_id = NULL
)

```

Arguments

subject	character of project and subject, such as "demo/YAB", or RAVESubject instance
electrodes	integer vector of electrodes, or a character that can be parsed by parse_svec
reference_name	reference name to be loaded
...	ignored
quiet	whether to quietly load the data
repository_id	used internally
epoch_name	epoch name to be loaded, or a RAVEEpoch instance
time_windows	a list of time windows that are relative to epoch onset time; need to pass the validation validate_time_window
stitch_events	events to stitch, default is NULL, meaning when loading data, the time is relative to the trial onset (column "Time" in the epoch file); set to a character of length 2, representing the events if time is not relative to trial onset. Please remove the prefix. For example, for a column named "Event_name", the event name is "name".
env	environment to evaluate
blocks	one or more session blocks to load
raw	whether to load from original (before processing) data; if true, then time-frequency data will not be loaded.
signal_type	electrode signal type (length of one) to be considered; default is 'LFP'. This option rarely needs to change unless you really want to check the power data from other types. For other signal types, check SIGNAL_TYPES
time_frequency	whether to load time-frequency data when preparing block data
verbose	whether to show progress

Value

A [fastmap2](#) (basically a list) of objects. Depending on the functions called, the following items may exist in the list:

```

subject A RAVESubject instance
epoch_name Same as input epoch_name
epoch A RAVEEpoch instance
reference_name Same as input reference_name
reference_table A data frame of reference
electrode_table A data frame of electrode information

```

frequency A vector of frequencies
 time_points A vector of time-points
 power_list A list of power data of the electrodes
 power_dimnames A list of trial indices, frequencies, time points, and electrodes that are loaded

progress_with_logger *Enhanced progress with logger message*

Description

For best performance, please install 'ravedash'. This function can replace [progress2](#).

Usage

```
progress_with_logger(
  title,
  max = 1,
  ...,
  quiet = FALSE,
  session = shiny::getDefaultReactiveDomain(),
  shiny_auto_close = FALSE,
  outputId = NULL,
  log
)
```

Arguments

title, max, ..., quiet, session, shiny_auto_close
 see [progress2](#)

outputId will be used if package 'shidashi' is installed, otherwise will be ignored

log function, NULL, or missing; default is missing, which will use logger function in the package 'ravedash', or [cat2](#) if 'ravedash' is not installed. If log=NULL, then the message will be suppressed in 'shiny' applications. If a function provided, then the function will be called.

Value

A list, see [progress2](#)

py_nipy_coreg

*Register 'CT' to 'MR' images via 'nipy' script***Description**

Align 'CT' using `nipy.algorithms.registration.histogram_registration`.

Usage

```
py_nipy_coreg(
    ct_path,
    mri_path,
    clean_source = TRUE,
    inverse_target = TRUE,
    precenter_source = TRUE,
    smooth = 0,
    reg_type = c("rigid", "affine"),
    interp = c("pv", "tri"),
    similarity = c("crl1", "cc", "cr", "mi", "nmi", "slr"),
    optimizer = c("powell", "steepest", "cg", "bfgs", "simplex"),
    tol = 1e-04,
    dry_run = FALSE
)

cmd_run_nipy_coreg(
    subject,
    ct_path,
    mri_path,
    clean_source = TRUE,
    inverse_target = TRUE,
    precenter_source = TRUE,
    reg_type = c("rigid", "affine"),
    interp = c("pv", "tri"),
    similarity = c("crl1", "cc", "cr", "mi", "nmi", "slr"),
    optimizer = c("powell", "steepest", "cg", "bfgs", "simplex"),
    dry_run = FALSE,
    verbose = FALSE
)
```

Arguments

<code>ct_path, mri_path</code>	absolute paths to 'CT' and 'MR' image files
<code>clean_source</code>	whether to replace negative 'CT' values with zeros; default is true
<code>inverse_target</code>	whether to inverse 'MRI' color intensity; default is true

precenter_source	whether to adjust the 'CT' transform matrix before alignment, such that the origin of 'CT' is at the center of the volume; default is true. This option may avoid the case that 'CT' is too far-away from the 'MR' volume at the beginning of the optimization
smooth, interp, optimizer, tol	optimization parameters, see 'nipy' documentation for details.
reg_type	registration type, choices are 'rigid' or 'affine'
similarity	the cost function of the alignment; choices are 'cr11' ('L1' regularized correlation), 'cc' (correlation coefficient), 'cr' (correlation), 'mi' (mutual information), 'nmi' (normalized mutual information), 'slr' (likelihood ratio). In reality I personally find 'cr11' works best in most cases, though many tutorials suggest 'nmi'.
dry_run	whether to dry-run the script and to print out the command instead of executing the code; default is false
subject	'RAVE' subject
verbose	whether to verbose command; default is false

Value

Nothing is returned from the function. However, several files will be generated at the 'CT' path:

'ct_in_t1.nii' aligned 'CT' image; the image is also re-sampled into 'MRI' space
 'CT_IJK_to_MR_RAS.txt' transform matrix from volume 'IJK' space in the original 'CT' to the 'RAS' anatomical coordinate in 'MR' scanner
 'CT_RAS_to_MR_RAS.txt' transform matrix from scanner 'RAS' space in the original 'CT' to 'RAS' in 'MR' scanner space

rave-raw-validation	<i>Validate raw files in 'rave' directory</i>
---------------------	---

Description

Validate subjects and returns whether the subject can be imported into 'rave'

Usage

```
validate_raw_file(
    subject_code,
    blocks,
    electrodes,
    format,
    data_type = c("continuous"),
    ...
)

IMPORT_FORMATS
```

Arguments

subject_code	subject code, direct folder under 'rave' raw data path
blocks	block character, direct folder under subject folder. For raw files following 'BIDS' convention, see details
electrodes	electrodes to verify
format	integer or character. For characters, run names(IMPORT_FORMATS)
data_type	currently only support continuous type of signals
...	other parameters used if validating 'BIDS' format; see details.

Format

An object of class list of length 7.

Details

Six types of raw file structures are supported. They can be basically classified into two categories: 'rave' native raw structure and 'BIDS-iEEG' structure.

In 'rave' native structure, subject folders are stored within the root directory, which can be obtained via `raveio_getopt('raw_data_dir')`. Subject directory is the subject code. Inside of subject folder are block files. In 'rave', term 'block' is the combination of session, task, and run. Within each block, there should be 'iEEG' data files.

In 'BIDS-iEEG' format, the root directory can be obtained via `raveio_getopt('bids_data_dir')`. 'BIDS' root folder contains project folders. This is unlike 'rave' native raw data format. Subject folders are stored within the project directories. The subject folders start with 'sub-'. Within subject folder, there are session folders with prefix 'ses-'. Session folders are optional. 'iEEG' data is stored in 'ieeg' folder under the session/subject folder. 'ieeg' folder should contain at least

electrodes.tsv sub-<label>*_{electrodes}.tsv

'iEEG' description sub-<label>*_{task}-<label>_{run}-<index>_{ieeg}.json

'iEEG' data file sub-<label>*_{task}-<label>_{run}-<index>_{ieeg}.<ext>, in current 'rave', only extensions '.vhdr+.eeg/.dat' ('BrainVision') or 'EDF' (or plus) are supported.

When format is 'BIDS', project_name must be specified.

The following formats are supported:

'mat/.h5 file per electrode per block' 'rave' native raw format, each block folder contains multiple 'Matlab' or 'HDF5' files. Each file corresponds to a channel/electrode. File names should follow 'xxx001.mat' or 'xxx001.h5'. The numbers before the extension are channel numbers.

'Single .mat/.h5 file per block' 'rave' native raw format, each block folder contains **only** one 'Matlab' or 'HDF5' file. The file name can be arbitrary, but extension must be either '.mat' or '.h5'. Within the file there should be a matrix containing all the data. The short dimension of the matrix will be channels, and larger side of the dimension corresponds to the time points.

'Single EDF(+) file per block' 'rave' native raw format, each block folder contains **only** one '.edf' file.

'Single BrainVision file (.vhdr+.eeg, .vhdr+.dat) per block' 'rave' native raw format, each block folder contains **only** two files. The first file is header '.vhdr' file. It contains all meta information. The second is either '.eeg' or '.dat' file containing the body, i.e. signal entries.

'BIDS & EDF(+)' 'BIDS' format. The data file should have '.edf' extension

'BIDS & BrainVision (.vhdr+.eeg, .vhdr+.dat)' 'BIDS' format. The data file should have '.vhdr'+'.eeg/.dat' extensions

Value

logical true or false whether the directory is valid. Attributes containing error reasons or snapshot of the data. The attributes might be:

snapshot	description of data found if passing the validation
valid_run_names	For 'BIDS' format, valid session+task+run name if passing the validation
reason	named list where the names are the reason why validation fails and values are corresponding sessions or electrodes or both.

rave-server	<i>Install and configure 'RAVE' server as background service using shiny-server</i>
-------------	---

Description

Works on 'Linux' and 'Mac' only.

Usage

```
rave_server_install(
  url = "https://github.com/rstudio/shiny-server/archive/refs/tags/v1.5.18.987.zip"
)

rave_server_configure(
  ports = 17283,
  user = Sys.info()[["user"]],
  rave_version = c("1", "2")
)
```

Arguments

url	'URL' to shiny-server 'ZIP' file to download
ports	integer vectors or character, indicating the port numbers to host 'RAVE' instances a valid port must be within the range from 1024 to 65535.
user	user to run the service as; default is the login user
rave_version	internally used; might be deprecated in the future

Value

nothing

Examples

```
## Not run:

# OS-specific. Please install R package `rpy2` first

# Install rave-server
rave_server_install()

# Let port 17283-17290 to host RAVE instance
rave_server_configure(ports = "17283-17290")

## End(Not run)
```

rave_brain	<i>Load 'FreeSurfer' or 'AFNI/SUMA' brain from 'RAVE'</i>
------------	---

Description

Create 3D visualization of the brain and visualize with modern web browsers

Usage

```
rave_brain(
  subject,
  surfaces = "pial",
  overlays = "aparc.a2009s+aseg",
  annotations = "label/aparc.a2009s",
  ...,
  use_template_if_missing = FALSE,
  include_electrodes = TRUE
)
```

Arguments

subject	character, list, or RAVESubject instance; for list or other objects, make sure <code>subject\$subject_id</code> is a valid 'RAVE' subject 'ID'
surfaces	one or more brain surface types from "pial", "white", "smoothwm", "pial-outer-smoothed", etc.; check threeBrain
overlays	volumes to overlay; default is 'aparc.a2009s+aseg'
annotations	surface annotation or curvature data to load; default is 'label/aparc.a2009s', referring to the '*h.aparc.a2009s.annot' under the label folder.
...	ignored, reserved for legacy code

usetemplateifmissing

whether to use template brain when the subject brain files are missing. If set to true, then a template (usually 'N27') brain will be displayed as an alternative solution, and electrodes will be rendered according to their 'MNI305' coordinates, or 'VertexNumber' if given.

include_electrodes

whether to include electrode in the model; default is true

Value

A 'threeBrain' instance if brain is found or usetemplateifmissing is set to true; otherwise returns NULL

Examples

```
# Please make sure DemoSubject is correctly installed
# The subject is ~1GB from Github

if(interactive()){
  brain <- rave_brain("demo/DemoSubject")

  if( !is.null(brain) ) { brain$plot() }
}
```

rave_command_line_path

Find and execute external command-line tools

Description

Find and execute external command-line tools

Usage

```
normalize_commandline_path(
  path,
  type = c("dcm2niix", "freesurfer", "fsl", "afni", "others"),
  unset = NA
)

cmd_dcm2niix(error_on_missing = TRUE, unset = NA)

cmd_freesurfer_home(error_on_missing = TRUE, unset = NA)
```

```

cmd_fsl_home(error_on_missing = TRUE, unset = NA)

cmd_afni_home(error_on_missing = TRUE, unset = NA)

cmd_homebrew(error_on_missing = TRUE, unset = NA)

is_dry_run()

```

Arguments

path	path to normalize
type	type of command
unset	default to return if the command is not found
error_on_missing	whether to raise errors if command is missing

Value

Normalized path to the command, or unset if command is missing.

rave_directories	<i>Returns a list of 'RAVE' directories</i>
------------------	---

Description

This function is internally used and should not be called directly.

Usage

```

rave_directories(
  subject_code,
  project_name,
  blocks = NULL,
  .force_format = c("", "native", "BIDS")
)

```

Arguments

subject_code	'RAVE' subject code
project_name	'RAVE' project name
blocks	session or block names, optional
.force_format	format of the data, default is automatically detected.

Value

A list of directories

rave_export	<i>Export 'RAVE' data</i>
-------------	---------------------------

Description

Export portable data for custom analyses.

Usage

```
rave_export(x, path, ...)

## Default S3 method:
rave_export(x, path, format = c("rds", "yaml", "json"), ...)

## S3 method for class 'rave_prepare_subject_raw_voltage_with_epoch'
rave_export(x, path, zip = FALSE, ...)

## S3 method for class 'rave_prepare_subject_voltage_with_epoch'
rave_export(x, path, zip = FALSE, ...)

## S3 method for class 'rave_prepare_power'
rave_export(x, path, zip = FALSE, ...)
```

Arguments

x	R object or 'RAVE' repositories
path	path to save to
...	passed to other methods
format	export format
zip	whether to zip the files

Value

Exported data path

Examples

```
x <- "my data"
path <- tempfile()
rave_export(x, path)

readRDS(path)

## Not run:
# Needs demo subject
path <- tempfile()
x <- prepare_subject_power("demo/DemoSubject")
```

```

# Export power data to path
rave_export(x, path)

## End(Not run)

```

rave_import

Import data into 'rave' projects

Description

Import files with predefined structures. Supported file formats include 'Matlab', 'HDF5', 'EDF(+)', 'BrainVision' ('.eeg/.dat/.vhdr'). Supported file structures include 'rave' native structure and 'BIDS' (very limited) format. Please see <https://openwetware.org/wiki/RAVE:ravepreprocess> for tutorials.

Usage

```

rave_import(
  project_name,
  subject_code,
  blocks,
  electrodes,
  format,
  sample_rate,
  conversion = NA,
  data_type = "LFP",
  task_runs = NULL,
  add = FALSE,
  ...
)

```

Arguments

project_name	project name, for 'rave' native structure, this can be any character; for 'BIDS' format, this must be consistent with 'BIDS' project name. For subjects with multiple tasks, see Section "'RAVE' Project"
subject_code	subject code in character. For 'rave' native structure, this is a folder name under raw directory. For 'BIDS', this is subject label without "sub-" prefix
blocks	characters, for 'rave' native format, this is the folder names subject directory; for 'BIDS', this is session name with "ses-". Section "Block vs. Session" for different meaning of "blocks" in 'rave' and 'BIDS'
electrodes	integers electrode numbers
format	integer from 1 to 6, or character. For characters, you can get options by running <code>names(IMPORT_FORMATS)</code>

sample_rate	sample frequency, must be positive
conversion	physical unit conversion, choices are NA, V, mV, uV
data_type	electrode signal type; see SIGNAL_TYPES
task_runs	for 'BIDS' formats only, see Section "Block vs. Session"
add	whether to add electrodes. If set to true, then only new electrodes are allowed to be imported, blocks will be ignored and trying to import electrodes that have been imported will still result in error.
...	other parameters

Value

None

'RAVE' Project

A 'rave' project can be very flexible. A project can refer to a task, a research objective, or "arbitrarily" as long as you find common research interests among subjects. One subject can appear in multiple projects with different blocks, hence project_name should be objective-based. There is no concept of "project" in 'rave' raw directory. When importing data, you choose subset of blocks from subjects forming a project.

When importing 'BIDS' data into 'rave', project_name must be consistent with 'BIDS' project name as a compromise. Once imported, you may change the project folder name in imported rave data directory to other names. Because once raw traces are imported, 'rave' data will become self-contained and 'BIDS' data are no longer required for analysis. This naming inconsistency will also be ignored.

Block vs. Session

'rave' and 'BIDS' have different definitions for a "chunk" of signals. In 'rave', we use "block". it means combination of session (days), task, and run, i.e. a block of continuous signals captured. Raw data files are supposed to be stored in file hierarchy of <raw-root>/<subject_code>/<block>/<datafiles>. In 'BIDS', sessions, tasks, and runs are separated, and only session names are indicated under subject folder. Because some previous compatibility issues, argument 'block' refers to direct folder names under subject directories. This means when importing data from 'BIDS' format, block argument needs to be session names to comply with 'subject/block' structure, and there is an additional mandatory argument task_runs especially designed for 'BIDS' format.

For 'rave' native raw data format, block will be as-is once imported.

For 'BIDS' format, task_runs will be treated as blocks once imported.

File Formats

Following file structure. Here use project "demo" and subject "YAB" and block "008"), electrode 14 as an example.

format=1, or ".mat/.h5 file per electrode per block" folder <raw>/YAB/008 contains 'Matlab' or 'HDF5' files per electrode. Data file name should look like "xxx_14.mat"

format=2, **or** "Single .mat/.h5 file per block" <raw>/YAB/008 contains only one 'Matlab' or 'HDF5' file. Data within the file should be a 2-dimensional matrix, where the column 14 is signal recorded from electrode 14

format=3, **or** "Single EDF(+) file per block" <raw>/YAB/008 contains only one 'edf' file

format=4, **or** "Single BrainVision file (.vhdr+.eeg, .vhdr+.dat) per block" <raw>/YAB/008 contains only one 'vhdr' file, and the data file must be inferred from the header file

format=5, **or** "BIDS & EDF(+)" <bids>/demo/sub-YAB/ses-008/ must contains *_electrodes.tsv, each run must have channel file. The channel files and electrode file must be consistent in names.

Argument task_runs is mandatory, characters, combination of session, task name, and run number. For example, a task header file in BIDS with name 'sub-YAB_ses-008_task-visual_run-01_ieeg.edf' has task_runs name as '008-visual-01', where the first '008' refers to session, 'visual' is task name, and the second '01' is run number.

format=6, **or** "BIDS & BrainVision (.vhdr+.eeg, .vhdr+.dat)" Same as previous format "BIDS & EDF(+)", but data files have 'BrainVision' formats.

rave_subject_format_conversion

Compatibility support for 'RAVE' 1.0 format

Description

Convert 'RAVE' subject generated by 2.0 pipeline such that 1.0 modules can use the data. The subject must have valid electrodes. The data must be imported, with time-frequency transformed to pass the validation before converting.

Usage

```
rave_subject_format_conversion(subject, verbose = TRUE, ...)
```

Arguments

subject	'RAVE' subject characters, such as 'demo/YAB', or a subject instance generated from RAVESubject
verbose	whether to verbose the messages
...	ignored, reserved for future use

Value

Nothing

RAVEAbstarctElectrode *Abstract definition of electrode class in RAVE*

Description

This class is not intended for direct use. Please create new child classes and implement some key methods.

Value

If simplify is enabled, and only one block is loaded, then the result will be a vector (type="voltage") or a matrix (others), otherwise the result will be a named list where the names are the blocks.

Public fields

subject subject instance ([RAVESubject](#))

number integer stands for electrode number or reference ID

reference reference electrode, either NULL for no reference or an electrode instance inherits RAVEAbstarctElectrode

epoch a [RAVEepoch](#) instance

stitch_events events to stitch, when loading window is not default to trial onset; must be NULL or a character vector of length 2

Active bindings

type signal type of the electrode, such as 'LFP', 'Spike', and 'EKG'; default is 'Unknown'

power_enabled whether the electrode can be used in power analyses such as frequency, or frequency-time analyses; this usually requires transforming the electrode raw voltage signals using signal processing methods such as 'Fourier', 'wavelet', 'Hilbert', 'multi-taper', etc. If an electrode has power data, then it's power data can be loaded via [prepare_subject_power](#) method.

is_reference whether this instance is a reference electrode

location location type of the electrode, see [LOCATION_TYPES](#) for details

exists whether electrode exists in subject

preprocess_file path to preprocess 'HDF5' file

power_file path to power 'HDF5' file

phase_file path to phase 'HDF5' file

voltage_file path to voltage 'HDF5' file

reference_name reference electrode name

epoch_name current epoch name

cache_root run-time cache path; NA if epoch or trial intervals are missing

trial_intervals trial intervals relative to epoch onset

Methods

Public methods:

- [RAVEAbstarctElectrode\\$new\(\)](#)
- [RAVEAbstarctElectrode\\$set_reference\(\)](#)
- [RAVEAbstarctElectrode\\$set_epoch\(\)](#)
- [RAVEAbstarctElectrode\\$clear_cache\(\)](#)
- [RAVEAbstarctElectrode\\$clear_memory\(\)](#)
- [RAVEAbstarctElectrode\\$load_data\(\)](#)
- [RAVEAbstarctElectrode\\$load_blocks\(\)](#)
- [RAVEAbstarctElectrode\\$clone\(\)](#)

Method `new()`: constructor

Usage:

```
RAVEAbstarctElectrode$new(subject, number, quiet = FALSE)
```

Arguments:

`subject` character or [RAVESubject](#) instance
`number` current electrode number or reference ID
`quiet` reserved, whether to suppress warning messages

Method `set_reference()`: set reference for instance

Usage:

```
RAVEAbstarctElectrode$set_reference(reference)
```

Arguments:

`reference` NULL or [RAVEAbstarctElectrode](#) instance instance

Method `set_epoch()`: set epoch instance for the electrode

Usage:

```
RAVEAbstarctElectrode$set_epoch(epoch, stitch_events = NULL)
```

Arguments:

`epoch` characters or [RAVEEpoch](#) instance. For characters, make sure "epoch_<name>.csv" is in meta folder.
`stitch_events` events to stitch, default is NULL, meaning when loading data, the time is relative to the trial onset (column "Time" in the epoch file); set to a character of length 2, representing the events if time is not relative to trial onset. Please remove the prefix. For example, for a column named "Event_name", the event name is "name".

Method `clear_cache()`: method to clear cache on hard drive

Usage:

```
RAVEAbstarctElectrode$clear_cache(...)
```

Arguments:

... implemented by child instances

Method `clear_memory()`: method to clear memory

Usage:

```
RAVEAbstarctElectrode$clear_memory(...)
```

Arguments:

... implemented by child instances

Method `load_data()`: method to load electrode data

Usage:

```
RAVEAbstarctElectrode$load_data(type)
```

Arguments:

type data type such as "power", "phase", "voltage", "wavelet-coefficient", or others depending on child class implementations

Method `load_blocks()`: load electrode block-wise data (with reference), useful when epoch is absent

Usage:

```
RAVEAbstarctElectrode$load_blocks(blocks, type, simplify = TRUE)
```

Arguments:

blocks session blocks

type data type such as "power", "phase", "voltage", "wavelet-coefficient".

simplify whether to simplify the result

Method `clone()`: The objects of this class are cloneable with this method.

Usage:

```
RAVEAbstarctElectrode$clone(deep = FALSE)
```

Arguments:

deep Whether to make a deep clone.

Examples

```
## Not run:

# To run this example, please download demo subject (~700 MB) from
# https://github.com/beauchamplab/rave/releases/tag/v0.1.9-beta

generator <- RAVEAbstarctElectrode

# load demo subject electrode 14
e <- generator$new("demo/DemoSubject", number = 14)

# set epoch
e$subject$epoch_names
e$set_epoch("auditory_onset")
head(e$epoch$table)

# set epoch range (-1 to 2 seconds relative to onset)
e$trial_intervals <- c(-1,2)
```

```

# or to set multiple ranges
e$trial_intervals <- list(c(-2,-1), c(0, 2))

# set reference
e$subject$reference_names
reference_table <- e$subject$meta_data(
  meta_type = "reference",
  meta_name = "default")
ref_name <- subset(reference_table, Electrode == 14)[["Reference"]]

# the reference is CAR type, mean of electrode 13-16,24
ref_name

# load & set reference
ref <- generator$new(e$subject, ref_name)
e$set_reference(ref)

## End(Not run)

```

RAVEEpoch

*Definition for epoch class***Description**

Trial epoch, contains the following information: Block experiment block/session string; Time trial onset within that block; Trial trial number; Condition trial condition. Other optional columns are Event_xxx (starts with "Event"). See <https://openwetware.org/wiki/RAVE:Epoching> or more details.

Value

```
self$table
```

If event is one of "trial onset", "default", "", or NULL, then the result will be "Time" column; if the event is found, then return will be the corresponding event column. When the event is not found and missing is "error", error will be raised; default is to return "Time" column, as it's trial onset and is mandatory.

If condition_type is one of "default", "", or NULL, then the result will be "Condition" column; if the condition type is found, then return will be the corresponding condition type column. When the condition type is not found and missing is "error", error will be raised; default is to return "Condition" column, as it's the default and is mandatory.

Public fields

name epoch name, character

subject RAVESubject instance

data a list of trial information, internally used

table trial epoch table

.columns epoch column names, internally used

Active bindings

columns columns of trial table
 n_trials total number of trials
 trials trial numbers
 available_events available events other than trial onset
 available_condition_type available condition type other than the default

Methods**Public methods:**

- [RAVEEpoch\\$new\(\)](#)
- [RAVEEpoch\\$trial_at\(\)](#)
- [RAVEEpoch\\$update_table\(\)](#)
- [RAVEEpoch\\$set_trial\(\)](#)
- [RAVEEpoch\\$get_event_colname\(\)](#)
- [RAVEEpoch\\$get_condition_colname\(\)](#)
- [RAVEEpoch\\$clone\(\)](#)

Method new(): constructor*Usage:*

```
RAVEEpoch$new(subject, name)
```

Arguments:

subject RAVESubject instance or character
 name character, make sure "epoch_<name>.csv" is in meta folder

Method trial_at(): get ith trial*Usage:*

```
RAVEEpoch$trial_at(i, df = TRUE)
```

Arguments:

i trial number
 df whether to return as data frame or a list

Method update_table(): manually update table field*Usage:*

```
RAVEEpoch$update_table()
```

Method set_trial(): set one trial*Usage:*

```
RAVEEpoch$set_trial(Block, Time, Trial, Condition, ...)
```

Arguments:

Block block string
 Time time in second

Trial positive integer, trial number
 Condition character, trial condition
 ... other key-value pairs corresponding to other optional columns

Method `get_event_colname()`: Get epoch column name that represents the desired event

Usage:

```
RAVEEpoch$get_event_colname(
  event = "",
  missing = c("warning", "error", "none")
)
```

Arguments:

event a character string of the event, see `$available_events` for all available events; set to "trial onset", "default", or blank to use the default
 missing what to do if event is missing; default is to warn

Method `get_condition_colname()`: Get condition column name that represents the desired condition type

Usage:

```
RAVEEpoch$get_condition_colname(
  condition_type,
  missing = c("warning", "error", "none")
)
```

Arguments:

condition_type a character string of the condition type, see `$available_condition_type` for all available condition types; set to "default" or blank to use the default
 missing what to do if condition type is missing; default is to warn if the condition column is not found.

Method `clone()`: The objects of this class are cloneable with this method.

Usage:

```
RAVEEpoch$clone(deep = FALSE)
```

Arguments:

deep Whether to make a deep clone.

Examples

```
# Please download DemoSubject ~700MB from
# https://github.com/beauchamplab/rave/releases/tag/v0.1.9-beta

## Not run:

# Load meta/epoch_auditory_onset.csv from subject demo/DemoSubject
epoch <-RAVEEpoch$new(subject = 'demo/DemoSubject',
  name = 'auditory_onset')

# first several trials
```

```

head(epoch$table)

# query specific trial
old_trial1 <- epoch$trial_at(1)

# Create new trial or change existing trial
epoch$set_trial(Block = '008', Time = 10,
                Trial = 1, Condition = 'AknownVmeant')
new_trial1 <- epoch$trial_at(1)

# Compare new and old trial 1
rbind(old_trial1, new_trial1)

# To get updated trial table, must update first
epoch$update_table()
head(epoch$table)

## End(Not run)

```

raveio-constants

The constant variables

Description

The constant variables

Usage

SIGNAL_TYPES

LOCATION_TYPES

MNI305_to_MNI152

Yael_IMAGE_TYPES

Format

An object of class character of length 5.

An object of class character of length 5.

An object of class matrix (inherits from array) with 4 rows and 4 columns.

An object of class character of length 10.

Details

SIGNAL_TYPES has the following options: 'LFP', 'Spike', 'EKG', 'Auxiliary', or 'Unknown'. As of 'raveio' 0.0.6, only 'LFP' (see [LFP_electrode](#)) signal type is supported.

LOCATION_TYPES is a list of the electrode location types: 'iEEG' (this includes the next two), 'sEEG' (stereo), 'ECoG' (surface), 'EEG' (scalp), 'Others'. See field 'location' in [RAVEAbstractElectrode](#)

MNI305_to_MNI152 is a 4-by-4 matrix converting 'MNI305' coordinates to 'MNI152' space. The difference of these two spaces is: 'MNI305' is an average of 305 human subjects, while 'MNI152' is the average of 152 people. These two coordinates differs slightly. While most of the 'MNI' coordinates reported by 'RAVE' and 'FreeSurfer' are in the 'MNI305' space, many other programs are expecting 'MNI152' coordinates.

RAVEMetaSubject	<i>Defines 'RAVE' subject class for meta analyses</i>
-----------------	---

Description

R6 class definition

Value

data frame

Super class

[raveio::RAVESubject](#) -> RAVEMetaSubject

Active bindings

project project instance of current subject; see [RAVEProject](#)
 project_name character string of project name
 subject_code character string of subject code
 subject_id subject ID: "project/subject"
 path subject root path
 rave_path 'rave' directory under subject root path
 meta_path meta data directory for current subject
 freesurfer_path 'FreeSurfer' directory for current subject. If no path exists, values will be NA
 preprocess_path preprocess directory under subject 'rave' path
 data_path data directory under subject 'rave' path
 cache_path path to 'FST' copies under subject 'data' path
 pipeline_path path to pipeline scripts under subject's folder
 note_path path that stores 'RAVE' related subject notes
 epoch_names possible epoch names

reference_names possible reference names
 reference_path reference path under 'rave' folder
 preprocess_settings preprocess instance; see [RAVEPreprocessSettings](#)
 blocks subject experiment blocks in current project
 electrodes all electrodes, no matter excluded or not
 raw_sample_rates voltage sample rate
 power_sample_rate power spectrum sample rate
 has_wavelet whether electrodes have wavelet transforms
 notch_filtered whether electrodes are Notch-filtered
 electrode_types electrode signal types

Methods

Public methods:

- [RAVEMetaSubject#print\(\)](#)
- [RAVEMetaSubject\\$new\(\)](#)
- [RAVEMetaSubject\\$meta_data\(\)](#)
- [RAVEMetaSubject\\$clone\(\)](#)

Method print(): override print method

Usage:

```
RAVEMetaSubject#print(...)
```

Arguments:

... ignored

Method new(): constructor

Usage:

```
RAVEMetaSubject$new(project_name, subject_code = NULL, strict = FALSE)
```

Arguments:

project_name character project name

subject_code character subject code

strict whether to check if subject folders exist

Method meta_data(): get subject meta data located in "meta/" folder

Usage:

```
RAVEMetaSubject$meta_data(  
  meta_type = c("electrodes", "frequencies", "time_points", "epoch", "references"),  
  meta_name = "default"  
)
```

Arguments:

meta_type choices are 'electrodes', 'frequencies', 'time_points', 'epoch', 'references'

```
meta_name if meta_type='epoch', read in 'epoch_<meta_name>.csv'; if meta_type='references',
  read in 'reference_<meta_name>.csv'.
```

Method `clone()`: The objects of this class are cloneable with this method.

Usage:

```
RAVEMetaSubject$clone(deep = FALSE)
```

Arguments:

`deep` Whether to make a deep clone.

See Also

[load_meta2](#)

RAVEPreprocessSettings

Defines preprocess configurations

Description

R6 class definition

Value

list of electrode type, number, etc.

NULL when no channel is composed. When `flat` is TRUE, a data frame of weights with the columns composing electrode channel numbers, composed channel number, and corresponding weights; if `flat` is FALSE, then a weight matrix;

Public fields

`current_version` current configuration setting version

`path` settings file path

`backup_path` alternative back up path for redundancy checks

`data` list of raw configurations, internally used only

`subject` [RAVESubject](#) instance

`read_only` whether the configuration should be read-only, not yet implemented

Active bindings

`version` configure version of currently stored files

`old_version` whether settings file is old format

`blocks` experiment blocks

`electrodes` electrode numbers

`sample_rates` voltage data sample rate

notch_filtered whether electrodes are notch filtered
 has_wavelet whether each electrode has wavelet transforms
 data_imported whether electrodes are imported
 data_locked whether electrode, blocks and sample rate are locked? usually when an electrode is imported into 'rave', that electrode is locked
 electrode_locked whether electrode is imported and locked
 electrode_composed composed electrode channels, not actual physically contacts, but is generated from those physically ones
 wavelet_params wavelet parameters
 notch_params Notch filter parameters
 electrode_types electrode signal types
 @freeze_blocks whether to free block, internally used
 @freeze_lfp_ecog whether to freeze electrodes that record 'LFP' signals, internally used
 @lfp_ecog_sample_rate 'LFP' sample rates, internally used
 all_blocks characters, all possible blocks even not included in some projects
 raw_path raw data path
 raw_path_type raw data path type, 'native' or 'bids'

Methods

Public methods:

- [RAVEPreprocessSettings\\$new\(\)](#)
- [RAVEPreprocessSettings\\$valid\(\)](#)
- [RAVEPreprocessSettings\\$has_raw\(\)](#)
- [RAVEPreprocessSettings\\$set_blocks\(\)](#)
- [RAVEPreprocessSettings\\$set_electrodes\(\)](#)
- [RAVEPreprocessSettings\\$set_sample_rates\(\)](#)
- [RAVEPreprocessSettings\\$migrate\(\)](#)
- [RAVEPreprocessSettings\\$electrode_info\(\)](#)
- [RAVEPreprocessSettings\\$save\(\)](#)
- [RAVEPreprocessSettings\\$get_compose_weights\(\)](#)

Method `new()`: constructor

Usage:

```
RAVEPreprocessSettings$new(subject, read_only = TRUE)
```

Arguments:

subject character or [RAVESubject](#) instance

read_only whether subject should be read-only (not yet implemented)

Method `valid()`: whether configuration is valid or not

Usage:

RAVEPreprocessSettings\$valid()

Method has_raw(): whether raw data folder exists

Usage:

RAVEPreprocessSettings\$has_raw()

Method set_blocks(): set blocks

Usage:

RAVEPreprocessSettings\$set_blocks(blocks, force = FALSE)

Arguments:

blocks character, combination of session task and run

force whether to ignore checking. Only used when data structure is not native, for example, 'BIDS' format

Method set_electrodes(): set electrodes

Usage:

```
RAVEPreprocessSettings$set_electrodes(
  electrodes,
  type = SIGNAL_TYPES,
  add = FALSE
)
```

Arguments:

electrodes integer vectors

type signal type of electrodes, see [SIGNAL_TYPES](#)

add whether to add to current settings

Method set_sample_rates(): set sample frequency

Usage:

RAVEPreprocessSettings\$set_sample_rates(srate, type = SIGNAL_TYPES)

Arguments:

srate sample rate, must be positive number

type electrode type to set sample rate. In 'rave', all electrodes with the same signal type must have the same sample rate.

Method migrate(): convert old format to new formats

Usage:

RAVEPreprocessSettings\$migrate(force = FALSE)

Arguments:

force whether to force migrate and save settings

Method electrode_info(): get electrode information

Usage:

RAVEPreprocessSettings\$electrode_info(electrode)

Arguments:

electrode integer

Method `save()`: save settings to hard disk

Usage:

`RAVEPreprocessSettings$save()`

Method `get_compose_weights()`: get weights of each composed channels

Usage:

`RAVEPreprocessSettings$get_compose_weights(flat = TRUE)`

Arguments:

`flat` whether to flatten the data frame; default is true

Examples

```
# The following example require downloading demo subject (~700 MB) from
# https://github.com/beauchamplab/rave/releases/tag/v0.1.9-beta

## Not run:

conf <- RAVEPreprocessSettings$new(subject = 'demo/DemoSubject')
conf$blocks # "008" "010" "011" "012"

conf$electrodes # 5 electrodes

# Electrode 14 information
conf$electrode_info(electrode = 14)

conf$data_imported # All 5 electrodes are imported

conf$data_locked # Whether block, sample rates should be locked

## End(Not run)
```

RAVEProject

Definition for 'RAVE' project class

Description

Definition for 'RAVE' project class

Definition for 'RAVE' project class

Value

character vector

true or false whether subject is in the project

A data table of pipeline time-stamps and directories

Active bindings

path project folder, absolute path

name project name, character

pipeline_path path to pipeline scripts under project's folder

Methods**Public methods:**

- `RAVEProject$print()`
- `RAVEProject$new()`
- `RAVEProject$subjects()`
- `RAVEProject$has_subject()`
- `RAVEProject$group_path()`
- `RAVEProject$subject_pipelines()`
- `RAVEProject$clone()`

Method `print()`: override print method

Usage:

`RAVEProject$print(...)`

Arguments:

... ignored

Method `new()`: constructor

Usage:

`RAVEProject$new(project_name, strict = TRUE)`

Arguments:

project_name character

strict whether to check project path

Method `subjects()`: get all imported subjects within project

Usage:

`RAVEProject$subjects()`

Method `has_subject()`: whether a specific subject exists in this project

Usage:

`RAVEProject$has_subject(subject_code)`

Arguments:

subject_code character, subject name

Method `group_path()`: get group data path for 'rave' module

Usage:

`RAVEProject$group_path(module_id, must_work = FALSE)`

Arguments:

module_id character, 'rave' module ID
 must_work whether the directory must exist; if not exists, should a new one be created?

Method subject_pipelines(): list saved pipelines

Usage:

```
RAVEProject$subject_pipelines(  
  pipeline_name,  
  cache = FALSE,  
  check = TRUE,  
  all = FALSE  
)
```

Arguments:

pipeline_name name of the pipeline
 cache whether to use cached registry
 check whether to check if the pipelines exist as directories
 all whether to list all pipelines; default is false; pipelines with the same label but older time-stamps will be hidden

Method clone(): The objects of this class are cloneable with this method.

Usage:

```
RAVEProject$clone(deep = FALSE)
```

Arguments:

deep Whether to make a deep clone.

RAVESubject

Defines 'RAVE' subject class

Description

R6 class definition

Value

data frame

integer vector of valid electrodes

The same as value

A named list of key-value pairs, or if one key is specified and simplify=TRUE, then only the value will be returned.

A data frame with four columns: 'namespace' for the group name of the entry (entries within the same namespace usually share same module), 'timestamp' for when the entry was registered. 'entry_name' is the name of the entry. If include_history is true, then multiple entries with the same 'entry_name' might appear since the obsolete entries are included. 'entry_value' is the value of the corresponding entry.

If `as_table` is `FALSE`, then returns as [RAVEEpoch](#) instance; otherwise returns epoch table; will raise errors when file is missing or the epoch is invalid.

If `simplify` is `true`, returns a vector of reference electrode names, otherwise returns the whole table; will raise errors when file is missing or the reference is invalid.

If `simplify` is `true`, returns a vector of electrodes that are valid (or won't be excluded) under given reference; otherwise returns a table. If `subset` is `true`, then the table will be subset and only rows with electrodes to be loaded will be kept.

If `simplify` is `true`, returns a vector of frequencies; otherwise returns a table.

A table of pipeline registry

A `PipelineTools` instance

Active bindings

`project` `project` instance of current subject; see [RAVEProject](#)

`project_name` character string of project name

`subject_code` character string of subject code

`subject_id` subject ID: "project/subject"

`path` subject root path

`rave_path` 'rave' directory under subject root path

`meta_path` meta data directory for current subject

`imaging_path` root path to imaging processing folder

`freesurfer_path` 'FreeSurfer' directory for current subject. If no path exists, values will be NA

`preprocess_path` preprocess directory under subject 'rave' path

`data_path` data directory under subject 'rave' path

`cache_path` path to 'FST' copies under subject 'data' path

`pipeline_path` path to pipeline scripts under subject's folder

`note_path` path that stores 'RAVE' related subject notes

`epoch_names` possible epoch names

`reference_names` possible reference names

`reference_path` reference path under 'rave' folder

`preprocess_settings` preprocess instance; see [RAVEPreprocessSettings](#)

`blocks` subject experiment blocks in current project

`electrodes` all electrodes, no matter excluded or not

`raw_sample_rates` voltage sample rate

`power_sample_rate` power spectrum sample rate

`has_wavelet` whether electrodes have wavelet transforms

`notch_filtered` whether electrodes are Notch-filtered

`electrode_types` electrode signal types

`electrode_composed` composed electrode channels, not actual physically contacts, but is generated from those physically ones

Methods

Public methods:

- `RAVESubject$print()`
- `RAVESubject$new()`
- `RAVESubject$meta_data()`
- `RAVESubject$valid_electrodes()`
- `RAVESubject$initialize_paths()`
- `RAVESubject$set_default()`
- `RAVESubject$get_default()`
- `RAVESubject$get_note_summary()`
- `RAVESubject$get_epoch()`
- `RAVESubject$get_reference()`
- `RAVESubject$get_electrode_table()`
- `RAVESubject$get_frequency()`
- `RAVESubject$list_pipelines()`
- `RAVESubject$load_pipeline()`
- `RAVESubject$clone()`

Method `print()`: override print method

Usage:

```
RAVESubject$print(...)
```

Arguments:

... ignored

Method `new()`: constructor

Usage:

```
RAVESubject$new(project_name, subject_code = NULL, strict = TRUE)
```

Arguments:

`project_name` character project name

`subject_code` character subject code

`strict` whether to check if subject folders exist

Method `meta_data()`: get subject meta data located in "meta/" folder

Usage:

```
RAVESubject$meta_data(
  meta_type = c("electrodes", "frequencies", "time_points", "epoch", "references"),
  meta_name = "default"
)
```

Arguments:

`meta_type` choices are 'electrodes', 'frequencies', 'time_points', 'epoch', 'references'

`meta_name` if `meta_type='epoch'`, read in 'epoch_<meta_name>.csv'; if `meta_type='references'`, read in 'reference_<meta_name>.csv'.

Method `valid_electrodes()`: get valid electrode numbers

Usage:

```
RAVESubject$valid_electrodes(reference_name, refresh = FALSE)
```

Arguments:

`reference_name` character, reference name, see `meta_name` in `self$meta_data` or [load_meta2](#) when `meta_type` is 'reference'

`refresh` whether to reload reference table before obtaining data, default is false

Method `initialize_paths()`: create subject's directories on hard disk

Usage:

```
RAVESubject$initialize_paths(include_freesurfer = TRUE)
```

Arguments:

`include_freesurfer` whether to create 'FreeSurfer' path

Method `set_default()`: set default key-value pair for the subject, used by 'RAVE' modules

Usage:

```
RAVESubject$set_default(key, value, namespace = "default")
```

Arguments:

`key` character

`value` value of the key

`namespace` file name of the note (without post-fix)

Method `get_default()`: get default key-value pairs for the subject, used by 'RAVE' modules

Usage:

```
RAVESubject$get_default(
  ...,
  default_if_missing = NULL,
  simplify = TRUE,
  namespace = "default"
)
```

Arguments:

... single key, or a vector of character keys

`default_if_missing` default value is any key is missing

`simplify` whether to simplify the results if there is only one key to fetch; default is TRUE

`namespace` file name of the note (without post-fix)

Method `get_note_summary()`: get summary table of all the key-value pairs used by 'RAVE' modules for the subject

Usage:

```
RAVESubject$get_note_summary(namespaces, include_history = FALSE)
```

Arguments:

`namespaces` namespaces for the entries; see method `get_default` or `set_default`. Default is all possible namespaces

include_history whether to include history entries; default is false

Method get_epoch(): check and get subject's epoch information

Usage:

```
RAVESubject$get_epoch(epoch_name, as_table = FALSE, trial_starts = 0)
```

Arguments:

epoch_name epoch name, depending on the subject's meta files

as_table whether to convert to [data.frame](#); default is false

trial_starts the start of the trial relative to epoch time; default is 0

Method get_reference(): check and get subject's reference information

Usage:

```
RAVESubject$get_reference(reference_name, simplify = FALSE)
```

Arguments:

reference_name reference name, depending on the subject's meta file settings

simplify whether to only return the reference column

Method get_electrode_table(): check and get subject's electrode table with electrodes that are load-able

Usage:

```
RAVESubject$get_electrode_table(
  electrodes,
  reference_name,
  subset = FALSE,
  simplify = FALSE
)
```

Arguments:

electrodes characters indicating integers such as "1-14,20-30", or integer vector of electrode numbers

reference_name see method get_reference

subset whether to subset the resulting data table

simplify whether to only return electrodes

Method get_frequency(): check and get subject's frequency table, time-frequency decomposition is needed.

Usage:

```
RAVESubject$get_frequency(simplify = TRUE)
```

Arguments:

simplify whether to simplify as vector

Method list_pipelines(): list saved pipelines

Usage:

```
RAVESubject$list_pipelines(
  pipeline_name,
  cache = FALSE,
  check = TRUE,
  all = FALSE
)
```

Arguments:

pipeline_name pipeline ID

cache whether to use cache registry to speed up

check whether to check if the pipelines exist

all whether to list all pipelines; default is false; pipelines with the same label but older time-stamps will be hidden

Method load_pipeline(): load saved pipeline

Usage:

```
RAVESubject$load_pipeline(directory)
```

Arguments:

directory pipeline directory name

Method clone(): The objects of this class are cloneable with this method.

Usage:

```
RAVESubject$clone(deep = FALSE)
```

Arguments:

deep Whether to make a deep clone.

See Also

[load_meta2](#)

read-brainvision-eeg *Load from 'BrainVision' file*

Description

Read in 'eeg' or 'ieeg' data from 'BrainVision' files with .eeg or .dat extensions.

Usage

```
read_eeg_header(file)
```

```
read_eeg_marker(file)
```

```
read_eeg_data(header, path = NULL)
```

Arguments

file	path to 'vhdr' header file
header	header object returned by read_eeg_header
path	optional, path to data file if original data file is missing or renamed; must be absolute path.

Details

A 'BrainVision' dataset is usually stored separately in header file (.vhdr), marker file (.vmrk, optional) and data file (.eeg or .dat). These files must store under a same folder to be read into R.

Header data contains channel information. Data "channel" contains channel name, reference, resolution and physical unit. "resolution" times digital data values is the physical value of the recorded data. read_eeg_data makes this conversion internally. "unit" is the physical unit of recordings. By default 'uV' means micro-volts.

Marker file that ends with .vmrk is optional. If the file is indicated by header file and exists, then a marker table will be included when reading headers. A marker table contains six columns: marker number, type, description, start position (in data point), size (duration in data points), and target channel (0 means applied for all channels).

Signal file name is usually contained within header file. Therefore it is desired that the signal file name never changed once created. However, in some cases when the signal files are renamed and cannot be indexed by header files, please specify path to force load signals from a different file.

Value

read_eeg_header returns a list containing information below:

raw	raw header contents
common	a list of descriptors of header
channels	table of channels, including number, reference, resolution and unit
sample_rate	sampling frequency
root_path	directory to where the data is stored
channel_counts	total channel counts
markers	NULL if marker file is missing, or list of marker description and table containing 6 columns.

read_eeg_data returns header, signal data and data description:

data	a matrix of signal values. Each row is a channel and each column is a time point.
------	---

Examples

```
header_file <- 'sub-01_ses-01_task-visual_run-01_ieeg.vhdr'

if( file.exists(header_file) ){
  # load a subject header
  header <- read_eeg_header(header_file)
```

```
# load entire signal
data <- read_eeg_data(header)

data$description
}
```

read-write-fst	<i>Read a 'fst' file</i>
----------------	--------------------------

Description

Read a 'fst' file

Usage

```
save_fst(x, path, ...)

load_fst(path, ..., as.data.table = TRUE)
```

Arguments

x	data frame to write to path
path	path to 'fst' file: must not be connection.
...	passed to read_fst or write_fst
as.data.table	passed to read_fst in fst package

read_csv_ieeg	<i>Read comma separated value file and ignore headers</i>
---------------	---

Description

Resolved some irregular 'iEEG' format where the header could be missing.

Usage

```
read_csv_ieeg(file, nrow = Inf, drop = NULL)
```

Arguments

file	comma separated value file to read from. The file must contains all numerical values
nrow	number of rows to read
drop	passed to fread

Details

The function checks the first two rows of comma separated value file. If the first row has different `storage.mode` than the second row, then the first row is considered header, otherwise header is treated missing. Note file must have at least two rows.

read_edf_header	<i>Read 'EDF(+)' or 'BDF(+)' file headers</i>
-----------------	---

Description

Wrapper of `readEdfHeader`, but added some information

Usage

```
read_edf_header(path)
```

Arguments

path	file path, passed to <code>readEdfHeader</code>
------	---

Details

The added names are: `isAnnot2`, `sampleRate2`, and `unit2`. To avoid conflict with other names, there is a "2" appended to each names. `isAnnot2` indicates whether each channel is annotation channel or recorded signals. `sampleRate2` is a vector of sample rates for each channels. `unit2` is physical unit of recorded signals. For 'iEEG' data, this is electric potential unit, and choices are 'V' for volt, 'mV' for millivolt, and 'uV' for micro-volt. For more details, see <https://www.edfplus.info/specs/edftexts.html>

Value

A list is header information of an 'EDF/BDF' file.

See Also

[readEdfHeader](#)

read_edf_signal	<i>Read 'EDF(+)' or 'BDF(+)' file signals</i>
-----------------	---

Description

Read 'EDF(+)' or 'BDF(+)' file signals

Usage

```
read_edf_signal(
    path,
    signal_numbers = NULL,
    convert_volt = c("NA", "V", "mV", "uV")
)
```

Arguments

path	file path, passed to readEdfHeader
signal_numbers	channel/electrode numbers
convert_volt	convert voltage (electric potential) to a new unit, NA means no conversion, other choices are 'V', 'mV', and 'uV'.

Value

A list containing header information, signal lists, and channel/electrode names. If signal_numbers is specified, the corresponding names should appear as selected_signal_names. get_signal() can get physical signals after unit conversion.

read_mat	<i>Read 'Matlab' files</i>
----------	----------------------------

Description

A compatible reader that can read both 'Matlab' files prior and after version 6.0

Usage

```
read_mat(file, ram = TRUE, engine = c("r", "py"))

read_mat2(file, ram = TRUE, engine = c("r", "py"))
```

Arguments

file	path to a 'Matlab' file
ram	whether to load data into memory. Only available when the file is in 'HDF5' format. Default is false and will load arrays, if set to true, then lazy-load data. This is useful when array is very large.
engine	method to read the file, choices are 'r' and 'py' ('Python'); if 'py' is chosen, make sure configure_conda is configured.

Details

[io_read_mat](#) can only read 'Matlab' files prior to version 6. After version 6, 'Matlab' uses 'HDF5' format to store its data, and `read_mat` can handle both cases.

The performance of `read_mat` can be limited when the file is too big or has many datasets as it reads all the data contained in 'Matlab' file into memory.

Value

A list of All the data stored in the file

See Also

[io_read_mat](#), [load_h5](#)

Examples

```
# Matlab .mat <= v7.3
x <- matrix(1:16, 4)
f <- tempfile()
ieegio::io_write_mat(list(x = x), con = f)

read_mat(f)

# Matlab .mat >= v7.3, using hdf5
# Make sure you have installed hdf5r
if( dipsaus::package_installed('hdf5r') ){

  f <- tempfile()
  save_h5(x, file = f, name = 'x')

  read_mat(f)

  # For v7.3, you don't have to load all data into RAM
  dat <- read_mat(f, ram = FALSE)
  dat

  dat$x[]

}
```

read_nsx_nev	<i>Read 'BlackRock' event and signal files</i>
--------------	--

Description

Current implementation supports minimum 2.3 file specification version. Please contact the package maintainer to add specification configurations if you want us to support older versions.

Usage

```
read_nsx_nev(
  paths,
  nev_path = NULL,
  header_only = FALSE,
  nev_data = TRUE,
  verbose = TRUE,
  ram = FALSE,
  force_update = FALSE,
  temp_path = file.path(tempdir(), "blackrock-temp")
)
```

Arguments

paths	'NSx' signal files, usually with file extensions such as '.ns1', '.ns2', '.ns3', '.ns4', '.ns5'.
nev_path	'NEV' event files, with file extension '.nev'
header_only	whether to load header information only and avoid reading signal arrays
nev_data	whether to load '.nev' comments and 'waveforms'
verbose	whether to print out progress when loading signal array
ram	whether to load signals into the memory rather than storing with filearray ; default is false
force_update	force updating the channel data even if the headers haven't changed
temp_path	temporary directory to store the channel data

`safe_read_csv`*Read comma separated value files with given column classes*

Description

Read comma separated value files with given column classes

Usage

```
safe_read_csv(  
  file,  
  header = TRUE,  
  sep = ",",  
  colClasses = NA,  
  skip = 0,  
  quote = "\"",  
  ...,  
  stringsAsFactors = FALSE  
)
```

Arguments

`file`, `header`, `sep`, `colClasses`, `skip`, `quote`, `stringsAsFactors`, ...
passed to `read.csv`

Details

Reading a comma separated value file using builtin function `read.csv` might result in some unexpected behavior. `safe_read_csv` does some preprocessing on the format so that it take cares of the following cases.

1. If `skip` exceeds the maximum rows of the data, return a blank data frame instead of raising error.
2. If row names are included in the file, `colClasses` automatically skip that column and starts from the second column
3. If length of `colClasses` does not equal to the number of columns, instead of cycling the class types, we set those columns to be NA type and let `read.csv` decide the default types.
4. `stringsAsFactors` is by default FALSE to be consistent with R 4.0, if the function is called in R 3.x.

Value

A data frame

Examples

```
f <- tempfile()
x <- data.frame(a = letters[1:10], b = 1:10, c = 2:11)

# ----- Auto-detect row names -----
# Write with rownames
utils::write.csv(x, f, row.names = LETTERS[2:11])

# read csv with base library utils
table1 <- utils::read.csv(f, colClasses = c('character', 'character'))

# 4 columns including row names
str(table1)

# read csv via safe_read_csv
table2 <- safe_read_csv(f, colClasses = c('character', 'character'))

# row names are automatically detected, hence 3 columns
# Only first columns are characters, the third column is auto
# detected as numeric
str(table2)

# read table without row names
utils::write.csv(x, f, row.names = FALSE)
table2 <- safe_read_csv(f, colClasses = c('character', 'character'))

# still 3 columns, and row names are 1:nrow
str(table2)

# ----- Blank data frame when nrow too large -----
# instead of raising errors, return blank data frame
safe_read_csv(f, skip = 1000)
```

safe_write_csv

*Save data to comma separated value files with backups***Description**

Save comma separated value files, if file exists, backup original file.

Usage

```
safe_write_csv(x, file, ..., quiet = FALSE)
```

Arguments

x, file, ...	pass to write.csv
quiet	whether to suppress overwrite message

Value

Normalized path of file

Examples

```
f <- tempfile()
x <- data.frame(a = 1:10)

# File not exists, same as write file, returns normalized `f`
safe_write_csv(x, f)

# Check whether file exists
file.exists(f)

# write again, and the old file will be copied
safe_write_csv(x, f)
```

save_h5

Save objects to 'HDF5' file without trivial checks

Description

Save objects to 'HDF5' file without trivial checks

Usage

```
save_h5(
  x,
  file,
  name,
  chunk = "auto",
  level = 4,
  replace = TRUE,
  new_file = FALSE,
  ctype = NULL,
  quiet = FALSE,
  ...
)
```

Arguments

x	an array, a matrix, or a vector
file	path to 'HDF5' file
name	path/name of the data; for example, "group/data_name"
chunk	chunk size
level	compress level from 0 - no compression to 10 - max compression

replace	should data be replaced if exists
new_file	should removing the file if old one exists
ctype	data type such as "character", "integer", or "numeric". If set to NULL then automatically detect types. Note for complex data please store separately the real and imaginary parts.
quiet	whether to suppress messages, default is false
...	passed to other LazyH5\$save

Value

Absolute path of the file saved

See Also

[load_h5](#)

Examples

```
file <- tempfile()
x <- array(1:120, dim = 2:5)

# save x to file with name /group/dataset/1
save_h5(x, file, '/group/dataset/1', chunk = dim(x))

# read data
y <- load_h5(file, '/group/dataset/1')
y[]
```

save_json	<i>Save or load R object in 'JSON' format</i>
-----------	---

Description

Save or load R object in 'JSON' format

Usage

```
save_json(
  x,
  con = stdout(),
  ...,
  digits = ceiling(-log10(.Machine$double.eps)),
  pretty = TRUE,
  serialize = TRUE
)

load_json(con, ..., map = NULL)
```

Arguments

x	R object to save
con	file or connection
...	other parameters to pass into toJSON or fromJSON
digits	number of digits to save
pretty	whether the output should be pretty
serialize	whether to save a serialized version of x; see 'Examples'.
map	a map to save the results

Value

save_json returns nothing; load_json returns an R object.

Examples

```
# Serialize
save_json(list(a = 1, b = function(){}))

# use toJSON
save_json(list(a = 1, b = function(){}), serialize = FALSE)

# Demo of using serializer
f1 <- tempfile(fileext = ".json")
save_json(x ~ y + 1, f1)

load_json(f1)

unlink(f1)
```

save_meta2

Function to save meta data to 'RAVE' subject

Description

Function to save meta data to 'RAVE' subject

Usage

```
save_meta2(data, meta_type, project_name, subject_code)
```

Arguments

data	data table
meta_type	see load meta
project_name	project name
subject_code	subject code

Value

Either none if no meta matched or the absolute path of file saved.

save_yaml	<i>Write named list to file</i>
-----------	---------------------------------

Description

Write named list to file

Usage

```
save_yaml(x, file, ..., sorted = FALSE)
```

Arguments

x	a named list, fastmap2 , or anything that can be transformed into named list via <code>as.list</code>
file, ...	passed to write_yaml
sorted	whether to sort the results by name; default is false

Value

Normalized file path

See Also

[fastmap2](#), [load_yaml](#), [read_yaml](#), [write_yaml](#)

Examples

```
x <- list(a = 1, b = 2)
f <- tempfile()

save_yaml(x, f)

load_yaml(f)

map <- dipsaus::fastmap2(missing_default = NA)
```

```
map$c <- 'lol'
load_yaml(f, map = map)

map$a
map$d
```

Spike_electrode	<i>Class definition for micro-wire spike channels</i>
-----------------	---

Description

Class definition for micro-wire spike channels

Class definition for micro-wire spike channels

Value

If simplify is enabled, and only one block is loaded, then the result will be a vector (type="voltage") or a matrix (others), otherwise the result will be a named list where the names are the blocks.

Super class

`raveio::RAVEAbstractElectrode` -> Spike_electrode

Active bindings

h5_fname 'HDF5' file name

valid whether current electrode is valid: subject exists and contains current electrode or reference;
subject electrode type matches with current electrode type

raw_sample_rate voltage sample rate

power_sample_rate power/phase sample rate

preprocess_info preprocess information

voltage_file path to voltage 'HDF5' file

Methods

Public methods:

- `Spike_electrode$print()`
- `Spike_electrode$set_reference()`
- `Spike_electrode$new()`
- `Spike_electrode$.load_noref_voltage()`
- `Spike_electrode$.load_raw_voltage()`
- `Spike_electrode$load_data()`
- `Spike_electrode$load_blocks()`

- `Spike_electrode$clear_cache()`
- `Spike_electrode$clear_memory()`
- `Spike_electrode$clone()`

Method `print()`: print electrode summary

Usage:

`Spike_electrode$print()`

Method `set_reference()`: set reference for current electrode

Usage:

`Spike_electrode$set_reference(reference)`

Arguments:

reference either NULL or LFP_electrode instance

Method `new()`: constructor

Usage:

`Spike_electrode$new(subject, number, quiet = FALSE)`

Arguments:

subject, number, quiet see constructor in [RAVEAbstarctElectrode](#)

Method `.load_noref_voltage()`: load non-referenced voltage (internally used)

Usage:

`Spike_electrode$.load_noref_voltage(reload = FALSE)`

Arguments:

reload whether to reload cache

srate voltage signal sample rate

Method `.load_raw_voltage()`: load raw voltage (no process)

Usage:

`Spike_electrode$.load_raw_voltage(reload = FALSE)`

Arguments:

reload whether to reload cache

Method `load_data()`: method to load electrode data

Usage:

`Spike_electrode$load_data(type = c("raw-voltage", "voltage"))`

Arguments:

type data type such as "power", "phase", "voltage", "wavelet-coefficient", and "raw-voltage".

For "power", "phase", and "wavelet-coefficient", 'Wavelet' transforms are required.

For "voltage", 'Notch' filters must be applied. All these types except for "raw-voltage" will be referenced. For "raw-voltage", no reference will be performed since the data will be the "raw" signal (no processing).

Method `load_blocks()`: load electrode block-wise data (with no reference), useful when epoch is absent

Usage:

```
Spike_electrode$load_blocks(
  blocks,
  type = c("raw-voltage", "voltage"),
  simplify = TRUE
)
```

Arguments:

`blocks` session blocks

`type` data type such as "power", "phase", "voltage", "raw-voltage" (with no filters applied, as-is from imported), "wavelet-coefficient". Note that if type is "raw-voltage", then the data only needs to be imported; for "voltage" data, 'Notch' filters must be applied; for all other types, 'Wavelet' transforms are required.

`simplify` whether to simplify the result

Method `clear_cache()`: method to clear cache on hard drive

Usage:

```
Spike_electrode$clear_cache(...)
```

Arguments:

`...` ignored

Method `clear_memory()`: method to clear memory

Usage:

```
Spike_electrode$clear_memory(...)
```

Arguments:

`...` ignored

Method `clone()`: The objects of this class are cloneable with this method.

Usage:

```
Spike_electrode$clone(deep = FALSE)
```

Arguments:

`deep` Whether to make a deep clone.

Tensor

R6 Class for large Tensor (Array) in Hybrid Mode

Description

can store on hard drive, and read slices of GB-level data in seconds

Value

self
the sliced data
a data frame with the dimension names as index columns and value_name as value column
original array
the collapsed data

Public fields

dim dimension of the array
dimnames dimension names of the array
use_index whether to use one dimension as index when storing data as multiple files
hybrid whether to allow data to be written to disk
last_used timestamp of the object was read
temporary whether to remove the files once garbage collected

Active bindings

varnames dimension names (read-only)
read_only whether to protect the swap files from being changed
swap_file file or files to save data to

Methods**Public methods:**

- `Tensor$do_finalize()`
- `Tensor$print()`
- `Tensor$.use_multi_files()`
- `Tensor$new()`
- `Tensor$subset()`
- `Tensor$flatten()`
- `Tensor$to_swap()`
- `Tensor$to_swap_now()`
- `Tensor$get_data()`
- `Tensor$set_data()`
- `Tensor$collapse()`
- `Tensor$operate()`

Method `do_finalize()`: release resource and remove files for temporary instances

Usage:

`Tensor$do_finalize()`

Method `print()`: print out the data dimensions and snapshot

Usage:

```
Tensor$print(...)
```

Arguments:

... ignored

Method `.use_multi_files()`: Internally used, whether to use multiple files to cache data instead of one

Usage:

```
Tensor$.use_multi_files(mult)
```

Arguments:

mult logical

Method `new()`: constructor

Usage:

```
Tensor$new(
  data,
  dim,
  dimnames,
  varnames,
  hybrid = FALSE,
  use_index = FALSE,
  swap_file = temp_tensor_file(),
  temporary = TRUE,
  multi_files = FALSE
)
```

Arguments:

data numeric array

dim dimension of the array

dimnames dimension names of the array

varnames characters, names of dimnames

hybrid whether to enable hybrid mode

use_index whether to use the last dimension for indexing

swap_file where to store the data in hybrid mode files to save data by index; default stores in `raveio_getopt('tensor_temp_path')`

temporary whether to remove temporary files when existing

multi_files if use_index is true, whether to use multiple

Method `subset()`: subset tensor

Usage:

```
Tensor$subset(..., drop = FALSE, data_only = FALSE, .env = parent.frame())
```

Arguments:

... dimension slices

drop whether to apply [drop](#) on subset data

data_only whether just return the data value, or wrap them as a Tensor instance

.env environment where ... is evaluated

Method `flatten()`: converts tensor (array) to a table (data frame)

Usage:

```
Tensor$flatten(include_index = FALSE, value_name = "value")
```

Arguments:

`include_index` logical, whether to include dimension names

`value_name` character, column name of the value

Method `to_swap()`: Serialize tensor to a file and store it via [write_fst](#)

Usage:

```
Tensor$to_swap(use_index = FALSE, delay = 0)
```

Arguments:

`use_index` whether to use one of the dimension as index for faster loading

`delay` if greater than 0, then check when last used, if not long ago, then do not swap to hard drive. If the difference of time is greater than delay in seconds, then swap immediately.

Method `to_swap_now()`: Serialize tensor to a file and store it via [write_fst](#) immediately

Usage:

```
Tensor$to_swap_now(use_index = FALSE)
```

Arguments:

`use_index` whether to use one of the dimension as index for faster loading

Method `get_data()`: restore data from hard drive to memory

Usage:

```
Tensor$get_data(drop = FALSE, gc_delay = 3)
```

Arguments:

`drop` whether to apply [drop](#) to the data

`gc_delay` seconds to delay the garbage collection

Method `set_data()`: set/replace data with given array

Usage:

```
Tensor$set_data(v)
```

Arguments:

`v` the value to replace the old one, must have the same dimension

notice the a tensor is an environment. If you change at one place, the data from all other places will change. So use it carefully.

Method `collapse()`: apply mean, sum, or median to collapse data

Usage:

```
Tensor$collapse(keep, method = "mean")
```

Arguments:

`keep` which dimensions to keep

method "mean", "sum", or "median"

Method `operate()`: apply the tensor by anything along given dimension

Usage:

```
Tensor$operate(
  by,
  fun = .Primitive("/"),
  match_dim,
  mem_optimize = FALSE,
  same_dimension = FALSE
)
```

Arguments:

by R object

fun function to apply

match_dim which dimensions to match with the data

mem_optimize optimize memory

same_dimension whether the return value has the same dimension as the original instance

Examples

```
if(!is_on_cran()){

# Create a tensor
ts <- Tensor$new(
  data = 1:18000000, c(3000,300,20),
  dimnames = list(A = 1:3000, B = 1:300, C = 1:20),
  varnames = c('A', 'B', 'C'))

# Size of tensor when in memory is usually large
# `lobstr::obj_size(ts)` -> 8.02 MB

# Enable hybrid mode
ts$to_swap_now()

# Hybrid mode, usually less than 1 MB
# `lobstr::obj_size(ts)` -> 814 kB

# Subset data
start1 <- Sys.time()
subset(ts, C ~ C < 10 & C > 5, A ~ A < 10)
#> Dimension: 9 x 300 x 4
#> - A: 1, 2, 3, 4, 5, 6,...
#> - B: 1, 2, 3, 4, 5, 6,...
#> - C: 6, 7, 8, 9
end1 <- Sys.time(); end1 - start1
#> Time difference of 0.188035 secs

# Join tensors
ts <- lapply(1:20, function(ii){
  Tensor$new(
```

```

    data = 1:9000, c(30,300,1),
    dimnames = list(A = 1:30, B = 1:300, C = ii),
    varnames = c('A', 'B', 'C'), use_index = 2)
  })
  ts <- join_tensors(ts, temporary = TRUE)

}

```

test_hdspeed	<i>Simple hard disk speed test</i>
--------------	------------------------------------

Description

Simple hard disk speed test

Usage

```

test_hdspeed(
  path = tempdir(),
  file_size = 1e+06,
  quiet = FALSE,
  abort_if_slow = TRUE,
  use_cache = FALSE
)

```

Arguments

path	an existing directory where to test speed, default is temporary local directory.
file_size	in bytes, default is 1 MB.
quiet	should verbose messages be suppressed?
abort_if_slow	abort test if hard drive is too slow. This usually happens when the hard drive is connected via slow internet: if the write speed is less than 0.1 MB per second.
use_cache	if hard drive speed was tested before, abort testing and return cached results or not; default is false.

Value

A vector of two: writing and reading speed in MB per seconds.

time_diff2	<i>Calculate time difference in seconds</i>
------------	---

Description

Calculate time difference in seconds

Usage

```
time_diff2(start, end, units = "secs", label = "")
```

Arguments

start, end	start and end of timer
units	passed to time_delta
label	rave-units label for display purpose.

Value

A number inherits rave-units class.

See Also

[as_rave_unit](#)

Examples

```
start <- Sys.time()
Sys.sleep(0.1)
end <- Sys.time()
dif <- time_diff2(start, end, label = 'Running ')
print(dif, digits = 4)

is.numeric(dif)

dif + 1
```

transform_point_to_template

Calculate template 'MNI' coordinates for points on native brain

Description

Calculate template 'MNI' coordinates for points on native brain

Usage

```
transform_point_to_template(
  subject,
  positions,
  space = c("scannerRAS", "tkrRAS"),
  mapping_method = c("volumetric", "surface"),
  flip_hemisphere = FALSE,
  verbose = TRUE,
  project_surface = "pial",
  volumetric_transform = c("auto", "affine", "nonlinear"),
  ...
)

transform_thinfilm_to_mni152(
  subject,
  flip_hemisphere = FALSE,
  interpolator = 0.3,
  n_segments = c(16, 16),
  group_labels = NULL,
  project_surface = "pial",
  volumetric_transform = c("auto", "affine", "nonlinear"),
  template_subject = c("cvs_avg35_inMNI152", "fsaverage", "bert", "MNI152")
)
```

Arguments

subject	'RAVE' subject
positions	optional matrix of 3 columns, either in scanner or surface space (specified by space); default is missing and will use the electrode localization results (electrodes.csv)
space	if positions is given, which native coordinate system should be used; default is native 'T1' (or 'scannerRAS'); alternative is 'FreeSurfer' surface coordinate (or 'tkrRAS')
mapping_method	whether the mapping is 'volumetric' or 'surface'; default is the former.
flip_hemisphere	whether to flip the hemisphere; default is FALSE
verbose	whether to verbose the mapping progress; default is true

project_surface	for surface mapping only, which surface to project electrodes onto; default is 'pial' surface, other common choices are 'white' for white-matter, or 'smoothwm' for smoothed white matter
volumetric_transform	for volume mapping only, which type of transform to use; default is 'auto' detecting and use non-linear deformation if exists, and fall back to 'affine' transform; other choices are 'affine' or 'nonlinear'
...	ignored
interpolator	whether the transform lean towards volume mapping (interpolator=0) or surface mapping (interpolator=1)
n_segments	positive integers with length of two: resolution of the mapping; default segments the thin-film array into 16 by 16 segments
group_labels	NULL (default) or a character vector indicating the group labels of thin-film electrodes; default assumes that all contacts are from thin-film electrodes.
template_subject	template subject to be mapped to; default is 'cvs_avg35_inMNI152', which is a 'MNI152' template generated by 'FreeSurfer'; other choices are 'fsaverage' and 'bert'

Value

A table of electrode 'MNI' coordinates.

Examples

```
if(interactive()) {  
  
  transform_point_to_template('demo/DemoSubject', mapping_method = "volumetric")  
  
}
```

url_neurosynth	<i>Get 'Neurosynth' website address using 'MNI152' coordinates</i>
----------------	--

Description

Get 'Neurosynth' website address using 'MNI152' coordinates

Usage

```
url_neurosynth(x, y, z)
```

Arguments

x, y, z	numerical values: the right-anterior-superior 'RAS' coordinates in 'MNI152' space
---------	---

Value

'Neurosynth' website address

validate_subject	<i>Validate subject data integrity</i>
------------------	--

Description

Check against existence, validity, and consistency

Arguments

subject	subject ID (character), or RAVESubject instance
method	validation method, choices are 'normal' (default) or 'basic' for fast checks; if set to 'normal', four additional validation parts will be tested (see parts with * in Section 'Value').
verbose	whether to print out the validation messages
version	data version, choices are 1 for 'RAVE' 1.0 data format, and 2 ('RAVE' 2.0 data format); default is 2

Value

A list of nested validation results. The validation process consists of the following parts in order:

Data paths (paths)

path the subject's root folder

path the subject's 'RAVE' folder (the 'rave' folder under the root directory)

raw_path the subject's raw data folder

data_path a directory storing all the voltage, power, phase data (before reference)

meta_path meta directory containing all the electrode coordinates, reference table, epoch information, etc.

reference_path a directory storing calculated reference signals

preprocess_path a directory storing all the preprocessing information

cache_path (**low priority**) data caching path

freesurfer_path (**low priority**) subject's 'FreeSurfer' directory

note_path (**low priority**) subject's notes

pipeline_path (**low priority**) a folder containing all saved pipelines for this subject

Preprocessing information (preprocess)

electrodes_set whether the subject has a non-empty electrode set

blocks_set whether the session block length is non-zero

sample_rate_set whether the raw sampling frequency is set to a valid, proper positive number

data_imported whether all the assigning electrodes have been imported
 notch_filtered whether all the 'LFP' and 'EKG' signals have been 'Notch' filtered
 has_wavelet whether all the 'LFP' signals are wavelet-transformed
 has_reference at least one reference has been generated in the meta folder
 has_epoch at least one epoch file has been generated in the meta folder
 has_electrode_file meta folder has electrodes.csv file

Meta information (meta)

meta_data_valid this item only exists when the previous preprocess validation is failed or incomplete
 meta_electrode_table the electrodes.csv file in the meta folder has correct format and consistent electrodes numbers to the preprocess information
 meta_reference_xxx (xxx will be replaced with actual reference names) checks whether the reference table contains all electrodes and whether each reference data exists
 meta_epoch_xxx (xxx will be replaced with actual epoch names) checks whether the epoch table has the correct formats and whether there are missing blocks indicated in the epoch files

Voltage data (voltage_data*)

voltage_preprocessing whether the raw preprocessing voltage data are valid. This includes data lengths are the same within the same blocks for each signal type
 voltage_data whether the voltage data (after 'Notch' filters) exist and readable. Besides, the lengths of the data must be consistent with the raw signals

Spectral power and phase (power_phase_data*)

power_data whether the power data exists for all 'LFP' signals. Besides, to pass the validation process, the frequency and time-point lengths must be consistent with the preprocess record
 power_data same as power_data but for the phase data

Epoch table (epoch_tables*) One or more sub-items depending on the number of epoch tables. To pass the validation, the event time for each session block must not exceed the actual signal duration. For example, if one session lasts for 200 seconds, it will invalidate the result if a trial onset time is later than 200 seconds.

Reference table (reference_tables*) One or more sub-items depending on the number of reference tables. To pass the validation, the reference data must be valid. The inconsistencies, for example, missing file, wrong frequency size, invalid time-point lengths will result in failure

validate_time_window *Validate time windows to be used*

Description

Make sure the time windows are valid intervals and returns a reshaped window list

Usage

```
validate_time_window(time_windows)
```

Arguments

time_windows vectors or a list of time intervals

Value

A list of time intervals (ordered, length of 2)

Examples

```
# Simple time window
validate_time_window(c(-1, 2))

# Multiple windows
validate_time_window(c(-1, 2, 3, 5))

# alternatively
validate_time_window(list(c(-1, 2), c(3, 5)))
validate_time_window(list(list(-1, 2), list(3, 5)))

## Not run:

# Incorrect usage (will raise errors)

# Invalid interval (length must be two for each intervals)
validate_time_window(list(c(-1, 2, 3, 5)))

# Time intervals must be in ascending order
validate_time_window(c(2, 1))

## End(Not run)
```

voltage_baseline	<i>Calculate voltage baseline</i>
------------------	-----------------------------------

Description

Calculate voltage baseline

Usage

```
voltage_baseline(
  x,
  baseline_windows,
  method = c("percentage", "zscore", "subtract_mean"),
```

```

    units = c("Trial", "Electrode"),
    ...
)

## S3 method for class 'rave_prepare_subject_raw_voltage_with_epoch'
voltage_baseline(
  x,
  baseline_windows,
  method = c("percentage", "zscore", "subtract_mean"),
  units = c("Trial", "Electrode"),
  electrodes,
  baseline_mean,
  baseline_sd,
  ...
)

## S3 method for class 'rave_prepare_subject_voltage_with_epoch'
voltage_baseline(
  x,
  baseline_windows,
  method = c("percentage", "zscore", "subtract_mean"),
  units = c("Trial", "Electrode"),
  electrodes,
  baseline_mean,
  baseline_sd,
  ...
)

## S3 method for class 'FileArray'
voltage_baseline(
  x,
  baseline_windows,
  method = c("percentage", "zscore", "subtract_mean"),
  units = c("Trial", "Electrode"),
  filebase = NULL,
  ...
)

## S3 method for class 'array'
voltage_baseline(
  x,
  baseline_windows,
  method = c("percentage", "zscore", "subtract_mean"),
  units = c("Trial", "Electrode"),
  ...
)

```

Arguments

x	R array, filearray , or 'rave_prepare_power' object created by prepare_subject_raw_voltage_with_epoch
baseline_windows	list of baseline window (intervals)
method	baseline method; choices are 'percentage' and 'zscore'; see 'Details' in baseline_array
units	the unit of the baseline; see 'Details'
...	passed to other methods
electrodes	the electrodes to be included in baseline calculation; for power repository object produced by prepare_subject_power only; default is all available electrodes in each of signal_types
baseline_mean, baseline_sd	internally used by 'RAVE' repository, provided baseline is not contained in the data. This is useful for calculating the baseline with data from other blocks.
filebase	where to store the output; default is NULL and is automatically determined

Details

The arrays must be three-mode tensor and must have valid named [dimnames](#). The dimension names must be 'Trial', 'Time', 'Electrode', case sensitive.

The `baseline_windows` determines the baseline windows that are used to calculate time-points of baseline to be included. This can be one or more intervals and must pass the validation function [validate_time_window](#).

The `units` determines the unit of the baseline. It can be either or both of 'Trial', 'Electrode'. The default value is both, i.e., baseline for each combination of trial and electrode.

Value

The same type as the inputs

Examples

```
## Not run:
# The following code need to download additional demo data
# Please see https://rave.wiki/ for more details

library(raveio)
repo <- prepare_subject_raw_voltage_with_epoch(
  subject = "demo/DemoSubject",
  time_windows = c(-1, 3),
  electrodes = c(14, 15))

##### Direct baseline on repository
voltage_baseline(
  x = repo, method = "zscore",
  baseline_windows = list(c(-1, 0), c(2, 3))
)
```

```

voltage_mean <- repo$raw_voltage$baseline$collapse(
  keep = c(1,3), method = "mean")
matplot(voltage_mean, type = "l", lty = 1,
  x = repo$raw_voltage$dimnames$Time,
  xlab = "Time (s)", ylab = "Voltage (z-scored)",
  main = "Mean coltage over trial (Baseline: -1~0 & 2~3)")
abline(v = 0, lty = 2, col = 'darkgreen')
text(x = 0, y = -0.5, "Aud-Onset ", col = "darkgreen", cex = 0.6, adj = c(1,1))

##### Alternatively, baseline on each electrode channel
voltage_mean2 <- sapply(repo$raw_voltage$data_list, function(inst) {
  re <- voltage_baseline(
    x = inst, method = "zscore",
    baseline_windows = list(c(-1, 0), c(2, 3)))
  rowMeans(re[])
})

# Same with floating difference
max(abs(voltage_mean - voltage_mean2)) < 1e-8

## End(Not run)

```

with_future_parallel *Enable parallel computing provided by 'future' package within the context*

Description

Enable parallel computing provided by 'future' package within the context

Usage

```

with_future_parallel(
  expr,
  env = parent.frame(),
  quoted = FALSE,
  on_failure = "multisession",
  max_workers = NA,
  ...
)

```

Arguments

expr	the expression to be evaluated
env	environment of the expr

quoted	whether expr has been quoted; default is false
on_failure	alternative 'future' plan to use if forking a process is disallowed; this usually occurs on 'Windows' machines; see details.
max_workers	maximum of workers; default is automatically set by <code>raveio_getopt("max_worker", 1L)</code>
...	additional parameters passing into make_forked_clusters

Details

Some 'RAVE' functions such as [prepare_subject_power](#) support parallel computing to speed up. However, the parallel computing is optional. You can enable it by wrapping the function calls within `with_future_parallel` (see examples).

The default plan is to use 'forked' R sessions. This is a convenient, fast, and relative simple way to create multiple R processes that share the same memories. However, on some machines such as 'Windows' the support has not yet been implemented. In such cases, the plan fall backs to a back-up specified by `on_failure`. By default, `on_failure` is 'multisession', a heavier implementation than forking the process, and slightly longer ramp-up time. However, the difference should be marginal for most of the functions.

When parallel computing is enabled, the number of parallel workers is specified by the option `raveio_getopt("max_worker", 1L)`.

Value

The evaluation results of expr

Examples

```
library(raveio)

demo_subject <- as_rave_subject("demo/DemoSubject", strict = FALSE)

if(dir.exists(demo_subject$path)) {
  with_future_parallel({
    prepare_subject_power("demo/DemoSubject")
  })
}
```

YAELOProcess

Class definition of 'YAELO' image pipeline

Description

Rigid-registration across multiple types of images, non-linear normalization from native brain to common templates, and map template atlas or 'ROI' back to native brain. See examples at [as_yael_process](#)

Value

whether the image has been set (or replaced)

Absolute path if the image

'RAVE' subject instance

Nothing

A list of moving and fixing images, with rigid transformations from different formats.

See method `get_template_mapping`

A list of input, output images, with forward and inverse transform files (usually two 'Affine' with one displacement field)

transformed image in 'ANTs' format

transformed image in 'ANTs' format

Nothing

A matrix of 3 columns, each row is a transformed points (invalid rows will be filled with NA)

A matrix of 3 columns, each row is a transformed points (invalid rows will be filled with NA)

Active bindings

`subject_code` 'RAVE' subject code

`image_types` allowed image types

`work_path` Working directory ('RAVE' imaging path)

Methods**Public methods:**

- `YAELProcess$new()`
- `YAELProcess$set_input_image()`
- `YAELProcess$get_input_image()`
- `YAELProcess$get_subject()`
- `YAELProcess$register_to_T1w()`
- `YAELProcess$get_native_mapping()`
- `YAELProcess$map_to_template()`
- `YAELProcess$get_template_mapping()`
- `YAELProcess$transform_image_from_template()`
- `YAELProcess$transform_image_to_template()`
- `YAELProcess$generate_atlas_from_template()`
- `YAELProcess$transform_points_to_template()`
- `YAELProcess$transform_points_from_template()`
- `YAELProcess$construct_ants_folder_from_template()`
- `YAELProcess$get_brain()`
- `YAELProcess$clone()`

Method `new()`: Constructor to instantiate the class

Usage:

```
YAELProcess$new(subject_code, image_types)
```

Arguments:

`subject_code` character code representing the subject

`image_types` vector of image types, such as 'T1w', 'CT', 'fGATIR'. All images except 'CT' will be considered 'preop' (before electrode implantation). Please use 'postop' to indicate if an image is taken after the implantation (for example, 'postopT1w')

Method `set_input_image()`: Set the raw input for different image types

Usage:

```
YAELProcess$set_input_image(
  path,
  type = YAEL_IMAGE_TYPES,
  overwrite = FALSE,
  on_error = c("warning", "error", "ignore")
)
```

Arguments:

`path` path to the image files in 'NIfTI' format

`type` type of the image

`overwrite` whether to overwrite existing images if the same type has been imported before; default is false

`on_error` when the file exists and `overwrite` is false, how should this error be reported; choices are 'warning' (default), 'error' (throw error and abort), or 'ignore'.

Method `get_input_image()`: Get image path

Usage:

```
YAELProcess$get_input_image(type = YAEL_IMAGE_TYPES)
```

Arguments:

`type` type of the image

Method `get_subject()`: Get 'RAVE' subject instance

Usage:

```
YAELProcess$get_subject(project_name = "YAEL", strict = FALSE)
```

Arguments:

`project_name` project name; default is 'YAEL'

`strict` passed to [as_rave_subject](#)

Method `register_to_T1w()`: Register other images to 'T1' weighted 'MRI'

Usage:

```
YAELProcess$register_to_T1w(image_type = "CT", reverse = FALSE, verbose = TRUE)
```

Arguments:

`image_type` type of the image to register, must be set via `process$set_input_image` first.

`reverse` whether to reverse the registration; default is false, meaning the fixed (reference) image is the 'T1'. When setting to true, then the 'T1' 'MRI' will become the moving image

`verbose` whether to print out the process; default is true

Method `get_native_mapping()`: Get the mapping configurations used by `register_to_T1w`

Usage:

```
YAELProcess$get_native_mapping(image_type = YAEL_IMAGE_TYPES, relative = FALSE)
```

Arguments:

`image_type` type of the image registered to 'T1' weighted 'MRI'

`relative` whether to use relative path (to the `work_path` field)

Method `map_to_template()`: Normalize native brain to 'MNI152' template

Usage:

```
YAELProcess$map_to_template(
  template_name = rpyants_builtin_templates(),
  use_images = c("T1w", "T2w", "T1wContrast", "fGATIR", "preopCT"),
  native_type = "T1w",
  verbose = TRUE,
  ...
)
```

Arguments:

`template_name` which template to use, choices are 'mni_icbm152_nlin_asym_09a', 'mni_icbm152_nlin_asym_09b', 'mni_icbm152_nlin_asym_09c', and 'fsaverage'.

`use_images` a vector of image types to use for normalization; default types are 'T1w', 'T2w', 'T1wContrast', 'fGATIR', and 'preopCT'. To use all available images for normalization, use wildcard "all"

`native_type` which type of image should be used to map to template; default is 'T1w'

`verbose` whether to print out the process; default is true

... additional tuning parameters passed to internal 'Python' code.

Method `get_template_mapping()`: Get configurations used for normalization

Usage:

```
YAELProcess$get_template_mapping(
  template_name = rpyants_builtin_templates(),
  native_type = "T1w",
  relative = FALSE
)
```

Arguments:

`template_name` which template is used

`native_type` which native image is mapped to template

`relative` whether the paths should be relative or absolute; default is false (absolute paths)

Method `transform_image_from_template()`: Apply transform from images (usually an atlas or 'ROI') on template to native space

Usage:

```
YAELOProcess$transform_image_from_template(
  template_roi_path,
  template_name = rpyants_builtin_templates(),
  native_type = "T1w",
  interpolator = c("auto", "nearestNeighbor", "linear", "gaussian", "bSpline",
    "cosineWindowedSinc", "welchWindowedSinc", "hammingWindowedSinc",
    "lanczosWindowedSinc", "genericLabel"),
  verbose = TRUE
)
```

Arguments:

`template_roi_path` path to the template image file which will be transformed into individuals' image

`template_name` templates to use

`native_type` which type of native image to use for calculating the coordinates (default 'T1w')

`interpolator` how to interpolate the 'voxels'; default is "auto": 'linear' for probabilistic map and 'nearestNeighbor' otherwise.

`verbose` whether the print out the progress

Method `transform_image_to_template()`: Apply transform to images (usually an atlas or 'ROI') from native space to template

Usage:

```
YAELOProcess$transform_image_to_template(
  native_roi_path,
  template_name = rpyants_builtin_templates(),
  native_type = "T1w",
  interpolator = c("auto", "nearestNeighbor", "linear", "gaussian", "bSpline",
    "cosineWindowedSinc", "welchWindowedSinc", "hammingWindowedSinc",
    "lanczosWindowedSinc", "genericLabel"),
  verbose = TRUE
)
```

Arguments:

`native_roi_path` path to the native image file that will be transformed into template

`template_name` templates to use

`native_type` which type of native image to use for calculating the coordinates (default 'T1w')

`interpolator` how to interpolate the 'voxels'; default is "auto": 'linear' for probabilistic map and 'nearestNeighbor' otherwise.

`verbose` whether the print out the progress

Method `generate_atlas_from_template()`: Generate atlas maps from template and morph to native brain

Usage:

```
YAELOProcess$generate_atlas_from_template(
  template_name = rpyants_builtin_templates(),
  atlas_folder = NULL,
  surfaces = NA,
  verbose = TRUE
)
```

Arguments:

template_name which template to use

atlas_folder path to the atlas folder (that contains the atlas files)

surfaces whether to generate surfaces (triangle mesh); default is NA (generate if not existed).

Other choices are TRUE for always generating and overwriting surface files, or FALSE to disable this function. The generated surfaces will stay in native 'T1' space.

verbose whether the print out the progress

Method transform_points_to_template(): Transform points from native images to template

Usage:

```
YAELProcess$transform_points_to_template(
  native_ras,
  template_name = rpyants_builtin_templates(),
  native_type = "T1w",
  verbose = TRUE
)
```

Arguments:

native_ras matrix or data frame with 3 columns indicating points sitting on native images in right-anterior-superior ('RAS') coordinate system.

template_name template to use for mapping

native_type native image type where the points sit on

verbose whether the print out the progress

Method transform_points_from_template(): Transform points from template images to native

Usage:

```
YAELProcess$transform_points_from_template(
  template_ras,
  template_name = rpyants_builtin_templates(),
  native_type = "T1w",
  verbose = TRUE
)
```

Arguments:

template_ras matrix or data frame with 3 columns indicating points sitting on template images in right-anterior-superior ('RAS') coordinate system.

template_name template to use for mapping

native_type native image type where the points sit on

verbose whether the print out the progress

Method construct_ants_folder_from_template(): Create a reconstruction folder (as an alternative option) that is generated from template brain to facilitate the '3D' viewer. Please make sure method map_to_template is called before using this method (or the program will fail)

Usage:

```
YAELOProcess$construct_ants_folder_from_template(
  template_name = rpyants_builtin_templates(),
  add_surfaces = TRUE
)
```

Arguments:

`template_name` template to use for mapping

`add_surfaces` whether to create surfaces that is morphed from template to local; default is TRUE. Please enable this option only if the cortical surfaces are not critical (for example, you are studying the deep brain structures). Always use 'FreeSurfer' if cortical information is used.

Method `get_brain()`: Get '3D' brain model

Usage:

```
YAELOProcess$get_brain(
  electrodes = TRUE,
  project_name = "YAELO",
  coord_sys = c("scannerRAS", "tkrRAS", "MNI152", "MNI305"),
  ...
)
```

Arguments:

`electrodes` whether to add electrodes to the viewers; can be logical, data frame, or a character (path to electrode table). When the value is TRUE, the electrode file under `project_name` will be loaded; when `electrodes` is a [data.frame](#), or path to a 'csv' file, then please specify `coord_sys` on what is the coordinate system used for columns "x", "y", and "z".

`project_name` project name under which the electrode table should be queried, if `electrodes=TRUE`

`coord_sys` coordinate system if `electrodes` is a data frame with columns "x", "y", and "z", available choices are 'scannerRAS' (defined by 'T1' weighted native 'MRI' image), 'tkrRAS' ('FreeSurfer' defined native 'TK-registered'), 'MNI152' (template 'MNI' coordinate system averaged over 152 subjects; this is the common "'MNI' coordinate space" we often refer to), and 'MNI305' (template 'MNI' coordinate system averaged over 305 subjects; this coordinate system used by templates such as 'fsaverage')

... passed to [threeBrain](#)

Method `clone()`: The objects of this class are cloneable with this method.

Usage:

```
YAELOProcess$clone(deep = FALSE)
```

Arguments:

`deep` Whether to make a deep clone.

Index

* datasets

- rave-raw-validation, [74](#)
- raveio-constants, [90](#)

- ants_coreg, [4](#)
- ants_morph_electrode (ants_coreg), [4](#)
- ants_mri_to_template (ants_coreg), [4](#)
- ants_preprocessing, [6](#)
- ants_registration, [5](#)
- archive_subject, [6](#)
- as_rave_project, [8](#)
- as_rave_subject, [9](#), [135](#)
- as_rave_unit, [9](#), [124](#)
- as_yael_process, [10](#), [133](#)
- Auxiliary_electrode, [11](#)

- backup_file, [13](#)
- baseline_array, [67](#), [131](#)
- BlackrockFile, [14](#)

- cache_path, [16](#)
- cache_root (cache_path), [16](#)
- cache_to_filearray, [17](#)
- cat2, [19](#), [72](#)
- catgl, [19](#)
- clear_cached_files (cache_path), [16](#)
- cmd-external (cmd_run_3dAllineate), [20](#)
- cmd_afni_home (rave_command_line_path), [78](#)
- cmd_dcm2niix (rave_command_line_path), [78](#)
- cmd_execute (cmd_run_3dAllineate), [20](#)
- cmd_freesurfer_home (rave_command_line_path), [78](#)
- cmd_fsl_home (rave_command_line_path), [78](#)
- cmd_homebrew (rave_command_line_path), [78](#)
- cmd_run_3dAllineate, [20](#)
- cmd_run_ants_coreg (ants_coreg), [4](#)
- cmd_run_ants_mri_to_template (ants_coreg), [4](#)
- cmd_run_dcm2niix (cmd_run_3dAllineate), [20](#)
- cmd_run_flirt (cmd_run_3dAllineate), [20](#)
- cmd_run_niftyreg_coreg (niftyreg_coreg), [64](#)
- cmd_run_nipy_coreg (py_nipy_coreg), [73](#)
- cmd_run_r (cmd_run_3dAllineate), [20](#)
- cmd_run_recon_all (cmd_run_3dAllineate), [20](#)
- cmd_run_recon_all_clinical (cmd_run_3dAllineate), [20](#)
- cmd_run_yael_preprocess, [23](#)
- collapse, [29](#)
- collapse2, [28](#)
- collapse_power, [26](#)
- compose_channel, [29](#)
- configure_conda, [108](#)
- convert_blackrock, [31](#)
- convert_electrode_table_to_bids, [32](#)

- data.frame, [102](#), [139](#)
- data.table, [35](#)
- dimnames, [67](#), [131](#)
- dir.create, [33](#)
- dir_create2, [33](#)
- download.file, [43](#)
- drop, [120](#), [121](#)

- ECoGTensor, [33](#), [67](#)
- export_table, [35](#)

- fastmap2, [62](#), [71](#), [115](#)
- filearray, [38](#), [67](#), [109](#), [131](#)
- filearray_create, [18](#)
- find_path, [36](#)
- findGlobals, [18](#)
- fread, [105](#)
- fromJSON, [114](#)

- generate_reference, [37](#)
- get_projects, [38](#)
- get_val2, [39](#)
- glue, [19](#)
- h5_names, [40](#)
- h5_valid, [40](#)
- import_electrode_table, [41](#)
- IMPORT_FORMATS (rave-raw-validation), [74](#)
- import_table (export_table), [35](#)
- ingest_regressor, [42](#)
- install_subject, [43](#)
- io_read_h5, [60](#)
- io_read_mat, [108](#)
- is.blank, [44](#)
- is.zerolenth, [44](#)
- is_dry_run (rave_command_line_path), [78](#)
- is_on_cran, [45](#)
- is_valid_ish, [39](#), [45](#)
- join_tensors, [47](#)
- lapply, [48](#)
- lapply_async, [48](#)
- lapply_async2, [48](#)
- LazyFST, [60](#)
- LazyH5, [60](#)
- LFP_electrode, [49](#), [91](#)
- LFP_reference, [54](#)
- load_bids_ieeg_header, [58](#)
- load_fst (read-write-fst), [105](#)
- load_fst_or_h5, [59](#)
- load_h5, [60](#), [108](#), [113](#)
- load_json (save_json), [113](#)
- load_meta2, [61](#), [93](#), [101](#), [103](#)
- load_yaml, [62](#), [115](#)
- LOCATION_TYPES, [84](#)
- LOCATION_TYPES (raveio-constants), [90](#)
- make_forked_clusters, [133](#)
- mgh_to_nii, [62](#)
- MNI305_to_MNI152 (raveio-constants), [90](#)
- mode, [46](#)
- new_electrode, [50](#), [63](#)
- new_reference, [38](#), [54](#)
- new_reference (new_electrode), [63](#)
- niftyreg_coreg, [64](#)
- normalize_commandline_path
(rave_command_line_path), [78](#)
- parse_svec, [63](#), [71](#)
- power_baseline, [66](#)
- prepare_subject_bare
(prepare_subject_bare0), [68](#)
- prepare_subject_bare0, [68](#)
- prepare_subject_phase
(prepare_subject_bare0), [68](#)
- prepare_subject_power, [42](#), [67](#), [84](#), [131](#),
[133](#)
- prepare_subject_power
(prepare_subject_bare0), [68](#)
- prepare_subject_raw_voltage_with_epoch,
[131](#)
- prepare_subject_raw_voltage_with_epoch
(prepare_subject_bare0), [68](#)
- prepare_subject_voltage_with_epoch
(prepare_subject_bare0), [68](#)
- prepare_subject_wavelet
(prepare_subject_bare0), [68](#)
- prepare_subject_with_blocks
(prepare_subject_bare0), [68](#)
- prepare_subject_with_epoch
(prepare_subject_bare0), [68](#)
- progress2, [72](#)
- progress_with_logger, [72](#)
- py_nipy_coreg, [73](#)
- rave-prepare (prepare_subject_bare0), [68](#)
- rave-raw-validation, [74](#)
- rave-server, [76](#)
- rave_brain, [77](#)
- rave_command_line_path, [78](#)
- rave_directories, [79](#)
- rave_export, [80](#)
- rave_import, [81](#)
- rave_server_configure (rave-server), [76](#)
- rave_server_install (rave-server), [76](#)
- rave_subject_format_conversion, [83](#)
- RAVEAbstarctElectrode, [12](#), [51](#), [55](#), [63](#), [84](#),
[91](#), [117](#)
- RAVEEpoch, [71](#), [84](#), [85](#), [87](#), [99](#)
- raveio-constants, [90](#)
- raveio::RAVEAbstarctElectrode, [11](#), [50](#),
[54](#), [116](#)
- raveio::RAVESubject, [91](#)
- raveio::Tensor, [34](#)

- raveio_getopt, [6](#), [48](#)
- RAVEMetaSubject, [91](#)
- RAVEPreprocessSettings, [92](#), [93](#), [99](#)
- RAVEProject, [8](#), [91](#), [96](#), [99](#)
- RAVESubject, [9](#), [10](#), [22](#), [38](#), [63](#), [71](#), [77](#), [83–85](#),
[93](#), [94](#), [98](#), [127](#)
- read-brainvision-eeg, [103](#)
- read-write-fst, [105](#)
- read.csv, [41](#)
- read_csv_ieeg, [105](#)
- read_edf_header, [106](#)
- read_edf_signal, [107](#)
- read_eeg_data (read-brainvision-eeg),
[103](#)
- read_eeg_header (read-brainvision-eeg),
[103](#)
- read_eeg_marker (read-brainvision-eeg),
[103](#)
- read_mat, [107](#)
- read_mat2 (read_mat), [107](#)
- read_nsx_nev, [109](#)
- read_yaml, [62](#), [115](#)
- readEdfHeader, [106](#)
- register_volume, [65](#)

- safe_read_csv, [110](#)
- safe_write_csv, [111](#)
- save_fst (read-write-fst), [105](#)
- save_h5, [61](#), [112](#)
- save_json, [113](#)
- save_meta2, [114](#)
- save_yaml, [62](#), [115](#)
- SIGNAL_TYPES, [30](#), [63](#), [71](#), [82](#), [95](#)
- SIGNAL_TYPES (raveio-constants), [90](#)
- Spike_electrode, [116](#)
- storage.mode, [106](#)
- system2, [22](#)

- Tensor, [29](#), [33](#), [34](#), [47](#), [118](#)
- test_hdspeed, [123](#)
- threeBrain, [77](#), [139](#)
- time_delta, [124](#)
- time_diff2, [124](#)
- toJSON, [114](#)
- transform_point_to_template, [125](#)
- transform_thinfilm_to_mni152
(transform_point_to_template),
[125](#)

- url_neurosynth, [126](#)

- validate_raw_file
(rave-raw-validation), [74](#)
- validate_subject, [127](#)
- validate_time_window, [67](#), [71](#), [128](#), [131](#)
- voltage_baseline, [129](#)

- with_future_parallel, [48](#), [132](#)
- write_fst, [105](#), [121](#)
- write_yaml, [62](#), [115](#)

- YAEI_IMAGE_TYPES (raveio-constants), [90](#)
- yael_preprocess, [6](#)
- yael_preprocess
(cmd_run_yael_preprocess), [23](#)
- YAEIProcess, [10](#), [133](#)

- zip, [7](#)